



COMBI 4000 M tympanografy & rhinomanometry







we are one of the leading companies in the field of ENT-diagnostics – worldwide –

COMBI 4000 M/Tymp + Rhino

The HOMOTH COMBI4000M PC module combines two ENT diagnostic and modern technology. Because of the double function of tympanografy and rhinomanometry you save costs and space.

In the Tymp-Mode, the HOMOTH Combi4000M diagnostic module it is able to do a full-automatic impedance measurement within a few seconds. Every function, as well as curves and data is presented at the PC monitor. The results can be verified before documentation to avoid misprintings.

The test automatic enables to operate the Tymp4000 with only one hand, so the head of the patient can be held steady with the other hand. For measurement the probe is held against the outer ear. A two coloured lamp on top of the handle indicates the correct fitting of the probe and the measurement is starting automatically. It is possible to measure the stapedius reflex ipsi- and contralateral. During this measurement, the middle ear pressure is held in the outer ear canal.

In the Rhino-Mode, time saving adaptation via nose plugs allows a direct use at the doctors working-place, so the measurement fits within the normal examination. All calculation of flow, pressure, percentage and resistance-coefficient is done in real time, so with completion of the measurement all results are present. Persons with beards, long hair and children with fear of masks, can be measured, too. For clinical use, a face-half-mask is available, too.

During the measurement the intensity of breathing can be controlled via two measuring bars. The measurement is carried out anterior under physiological conditions of self-breathing and allows a quantitatively objective statement of the resistance behaviour of the nose.

Applications:

- € Proof of allergies after provocation
 € Function-control after nose operation
- € Diagnostics at handicapped nose breathing
- € Control after dispensation of medicines

All calculation of flow, pressure, percentage and resistance coefficient are done in real time, so that with the end of the measurement the complete results are present. All functions, as well as curves and data are presented at the PC monitor. The results can be verified before documentation to avoid miss-printings. The data and measured curves can be stored at the PC or transmitted to a practice software.

Technical data rhinomanometry

measuring-methode: anterior self-breathing impedance and reflex measurement system: EN 60601 - 1 / 1-1 / 1-2 / MPG standards: EN 60601 - 1 / 1-1 / 1-2 / MPG standards: probetone: 226 Hz, 85 dB SPL airflow: 0 - 900 ml /s in- und exspiration 0 - 50 daPa difference pressure: + 200 to -400 daPa pressure range: via 2 measuring bars in the display function control: pressure delta: 300 daPa / sec 2 sec. for compliance averaging: via max 5 flow-curves measure time: pressure calibration: automatical before start 500, 1000, 2000 and 4000 Hz representation: as diagram and numeric reflex tones: at 85, 95 und 105 dB HL aut. calculation of: 1. pressure in daPa 2. flow in ml /sec sequence of tones: automatically reflex recognition: automatically 3. total-flow 4. percentage reflexes: ipsi and contralateral 5. resistance-coefficients pressure calibration: automatically at start 6. Point 2-5 at 75, 150 and 300 daPa state indication: 1. LED 3-colors in probe 2. detailed display symbols accessories: 12 nose-plugs (olives) 2 flow-probes changeover right/left: automatic or manual 2 sets of hoses mains 230 V 1 foot switch pc interface USB 1 probe holder air pump: very quiet syringe pump 1 instruction manual accessories: 1 probe with cable face-halfmask 1 set earplugs special accessories: 1 headphone DT 48 A 1 probe holder -----Pentium PC min. 500 MHz, system requirements: 1 instruction manual USB 2.0 port, Window 98, dimensions w= 290 / d= 270 / h= 125 mm ME, 2000, XP for technical modifications all rights reserved 0124

Rev. 04 / Stand 01-2007

HOMOTH Medizinelektronik GmbH & Co KG, Baumacker 1a, 22523 Hamburg Tel.: ..49..40-54 72 65 0 Fax: ..49..40-54 72 65 50 Internet //www.homoth.de e-mail: info@homoth.de

Technical data tympanografy