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Guide Wire

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### Soft tissue procedures

- Gingival Troughing for Crown Impressions
- Gingivectomy & Gingivoplasty
- Gingival Incision & Excision
- Soft-Tissue Crown Lengthening
- Hemostasis & Coagulation
- Excisional & Incisional Biopsies
- Exposure of Unerupted Teeth
- Fibromal Removal
- Frenectomy & Frenotomy
- Implant Recovery
- Incision & Drainage of Abscess
- Leukoplakia
- Pulpotomy as an Adjunct to Root Canal Therapy
- Operculectomy
- Oral Papillectomies
- Reduction of Gingival Hypertrophy
- Vestibuloplasty
- Treatment of Canker Sores, Herpetic & Aphthous Ulcers of the Oral Mucosa

### Peridontal procedures

- Sulcular Debridement (Removal of Diseased, Infected, Inflamed, & Necrosed Soft-Tissue in the Periodontal Pocket to Improve Clinical Indices Including Gingival Index, Gingival Bleeding Index, Probe Depth, Attachment Loss, & Tooth Mobility)
- Laser Soft-Tissue Curettage
- Laser Removal of Diseased, Infected, Inflamed & Necrosed Soft-Tissue Within the Periodontal Pocket
- Removal of Highly Inflamed Edematous Tissue Affected by Bacteria Penetration of the Pocket Lining & Junctional Epithelium

### Laser teeth whitening

- Laser Assisted Whitening/Bleaching of Teethure.

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Periodontal Therapy Gingivitis - 2W Periodontitis – Pulse, 500/500ms, 1-2.5W Operculectomy – CW, 3-3.5W Frenectomy – CW, 2-2.5W Gingivectomy – CW, 1.5-2W	Endodontics Root canal sterilization- CW, 1W Periapical abscess- CW 1.2-1.5W Caries-Pulse 1.2W, 500/500ms Desensitization- Pulse, 30/30ms, 3W Laser analgesia,-CW, 2W	Oral surgery Soft tissue excision-2W Ablation-2.5W Implant exposure-CW, 2-3W Peri-implantitis-2W Apthae and herpes- 1.5W
Aesthetic procedures Bleaching – CW, 1.5-3W Depigmentation – CW, 2-3.5W Gummy smile correction -2W Apthae and herpes, CW, 1.5-2W Crown lengthening, CW, 1.5W	Biostimulation Teeth desensitization –Pulse, 30/30ms, 0.5-5W Post extraction-0.8W TMJ -0.8W Haemostats-CW, 0.8-1.5W Leukoplakia-2W	

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tel. +48 74 843 81 09 fax +48 74 840 17 33

web: <https://www.videomed.eu> e-mail: [info@videomed.eu](mailto:info@videomed.eu)