

STERILE CONTAINER SYSTEMS BIO-BARRIER MODEL



1/1 FULL SIZE



3/4 QUARTER SIZE



1/2 HALF SIZE



EXTRA LARGE SIZE



EXTRA LONG SIZE

Aygün Bio-Barrier model sterilization container systems are designed with mechanical valves both in bottom and lid which run in a way that they do not need any disposable filter use and can be used minimum 5000 cycles. The mechanical valve systems set on the lid or box run with atmospheric pressure occurring in the autoclave during sterilization. When the sterilization process starts, because of the occurring (plus) pressure, the valves open and let the steam penetration into the sterilization container. When the vacuum process starts at the end of the sterilization, because of the occurring (minus) pressure, the steam outlet valve opens and the steam exhausts from the box. When the sterilization process is completed, as there is no pressure (plus or minus) in autoclave, both valves are in closed position. By this way, it gains qualification to keep the products in it as sterilized. The valves are ensured to be closed %100 due to the special silicone gaskets on mechanical valve systems. Moreover, the mechanical valve system these sterilization containers is demountable by the user in order to control it when desired.

Aygün Bio-Barrier model sterilization container systems are manufactured as five different length, width and height and can be supplied with safety lids for effective external protection. These safety lids do not interfere with sterilization and they provides protection to the valve mechanism against liquids and dust as a second lid. All models and sizes have specific baskets and accessories which can be viewed in following parts. Modular auto lock mechanism (thermo lock) promising sterility assurance and indicator window providing visibility of the indicator card result by visual control without opening the lid are also available for this sterilization container model.

STERILE CONTAINER SYSTEMS





Bio-Barrier Model
















































































» 3/4 Quarter Size

3/4 BIO-BARRIER MODEL

W: 285 mm

L: 475 mm




Specifications ➔	Lid (with barrier)	Lid (with barrier)	Safety Lid Lid (with barrier)	Safety Lid Lid (with barrier)
				
	Box (Non Perforated)	Box (with barrier)	Box (Non Perforated)	Box (with barrier)
	K7-240 	K7-260 	K7-350 	K7-370 
	K7-241 	K7-261 	K7-351 	K7-371 
	K7-242 	K7-262 	K7-352 	K7-372 
	K7-243 	K7-263 	K7-353 	K7-373 
	K7-244 	K7-264 	K7-354 	K7-374 
	K7-245 	K7-265 	K7-355 	K7-375 
	K7-246 	K7-266 	K7-356 	K7-376 
	K7-247 	K7-267 	K7-357 	K7-377 
	K7-248 	K7-268 	K7-358 	K7-378 
	K7-249 	K7-269 	K7-359 	K7-379 
	K7-250 	K7-270 	K7-360 	K7-380 
	K7-251 	K7-271 	K7-361 	K7-381 
	K7-252 	K7-272 	K7-362 	K7-382 
	K7-253 	K7-273 	K7-363 	K7-383 
	K7-254 	K7-274 	K7-364 	K7-384 
	K7-255 	K7-275 	K7-365 	K7-385 
	K7-256 	K7-276 	K7-366 	K7-386 
	K7-257 	K7-277 	K7-367 	K7-387 

➔ Ps: All container measurements cover the outer size of the lids.


BOX (with barrier)




outside dimension



outside dimension



inside dimension




code	outside dimension	inside dimension	weight
K7-300	475 x 272 h:100 mm	422 x 258 h:60 mm	2024 g
K7-301	475 x 272 h:116 mm	422 x 258 h:76 mm	2129 g
K7-302	475 x 272 h:138 mm	422 x 258 h:98 mm	2257 g
K7-303	475 x 272 h:185 mm	422 x 258 h:145 mm	2532 g


BOX (non perforated)




outside dimension



outside dimension



inside dimension



code	outside dimension	inside dimension	weight
K6-305	475 x 272 h:100 mm	422 x 258 h:94 mm	1609 g
K6-306	475 x 272 h:116 mm	422 x 258 h:110 mm	1714 g
K6-307	475 x 272 h:138 mm	422 x 258 h:132 mm	1842 g
K6-308	475 x 272 h:185 mm	422 x 258 h:179 mm	2117 g

LID-BARRIER SYSTEM

Dimension: 477 x 285 mm



code	color	weight
K7-340		1212 g
K7-341		
K7-342		
K7-343		
K7-344		
K7-345		

SAFETY LIDS

Dimension: 483 x 290 mm



code	color	weight
K7-406		998 g
K7-407		
K7-408		
K7-409		
K7-410		
K7-411		

STANDARD MODEL 29

BIO-BARRIER MODEL 63

Polymer's PPSU Lids 97

WINDOW MODEL 119

B₃ Barrier MODEL 139

PLASMA MODEL 173

A1 MODEL 185

Dental, Mini, Special 205


BASKETS 223

Autolock Mechanism Emulating Indicator 243

ACCESSORIES 251

BASKET MADE FROM STAINLESS STEEL WIRE MESH


WIRE BASKET



code	outside dimension
K2-1500	405 x 245 x h:30 mm
K2-1501	405 x 245 x h:50 mm
K2-1502	405 x 245 x h:70 mm
K2-1503	405 x 245 x h:100 mm



WIRE BASKET WITH LID



code	outside dimension
K3-1030	405 x 245 h:35 mm
K3-1031	405 x 245 h:55 mm
K3-1032	405 x 245 h:75 mm
K3-1033	405 x 245 h:105 mm
K3-1034	ONLY LID, for baskets




WIRE BASKET WITH PERFORATED METAL PLATE

WIRE BASKET



code	outside dimension
K2-1501 P	405 x 246 x h:50 mm
K2-1502 P	405 x 246 x h:70 mm
K2-1503 P	405 x 246 x h:100 mm



WIRE BASKET WITH LID



code	outside dimension
K3-1031 P	405 x 246 x h:55 mm
K3-1032 P	405 x 246 x h:75 mm
K3-1033 P	405 x 246 x h:105 mm
K3-1034 P	ONLY LID, for baskets






BASKET MADE FROM STAINLESS STEEL SHEET

BASKET (with stand)



	code	outside dimension
Box and Lid	K3-728	411 x 259 h:44 mm
	K3-729	411 x 259 h:64 mm
	K3-730	411 x 259 h:84 mm
	K3-731	411 x 259 h:114 mm
Box	K3-748	405 x 253 h:40 mm
	K3-749	405 x 253 h:60 mm
	K3-750	405 x 253 h:80 mm
	K3-751	405 x 253 h:110 mm
Lid	K3-762	ONLY LID, for baskets



The baskets which are without lid and less than 44 mm in height do not include the handles while the ones with lid include the handles.

STAND




code	height
K3-767	20 mm

STAND



code	height
K3-768	10 mm

STAND



code	height
K3-766	5 mm

IDENTIFICATION (aluminium)

LABEL HOLDER



code	color
K2-560	Red
K2-561	Blue
K2-562	Green
K2-563	Yellow
K2-564	Grey
K2-565	Black

LABELS



code	color
K2-570	Red
K2-571	Blue
K2-572	Green
K2-573	Yellow
K2-574	Grey
K2-575	Black

Coloured identification labels Engraved to your specification



Grey coloured recommended for use in Sterrad. Color fading By effect of Hydrogen Peroxide Gas Plasma on other colors.

STANDARD MODEL 29

BIO-BARRIER MODEL 63

Polymer's PPSU Lids 97

WINDOW MODEL 119

B₃ Barrier MODEL 139

PLASMA MODEL 173

A1 MODEL 185


Dental, Mini, Special 205

BASKETS 223

Autolock Mechanism Emulating Indicator 243

ACCESSORIES 251

DIVIDERS




K2-1540
405 x 250 x 80 mm
K2-1541
405 x 250 x 110 mm
K2-1542
405 x 250 x 130 mm

WIRE BASE



K2-1522
405 x 255 x 30 mm


DIVIDERS



K2-1530
405 x 250 x 80 mm
K2-1531
405 x 250 x 110 mm
K2-1532
405 x 250 x 130 mm

(-)

DRAPE RETAINER



K2-1521
405 x 255 x 30 mm

DRAPE RETAINER



K2-1520
405 x 255 mm

ALUMINIUM BASKET (with lid)



K3-1119
245 x 130 x 55 mm

ALUMINIUM BASKET (with lid)




K3-1120
245 x 130 x 70 mm

ALUMINIUM BASKET (with lid)




K3-1121
245 x 130 x 105 mm

SILICONE MAT



K2-524
380 x 230 mm
K2-525
440 x 230 mm

SILICONE MESH MAT



K2-534
395 x 245 mm

SILICONE HANDLE




K2-710
K2-711
K2-712
K2-713
K2-714
K2-715

IDENTIFICATION LABELS (aluminium)




K2-700
K2-701
K2-702
K2-703
K2-704
K2-705

LABELS (paper)



K2-673
100 pcs. With indicator.
For small size and special containers

LABELS (paper)



K2-672
100 pcs. With indicator.
For full and Half size containers


(-)

AUTO LOCK MECHANISM



K1-900
(2 Pcs)

EMULATING INDICATOR CARD (class 6)



K2-674
Paper 250 Pcs.



SECURITY SEAL (blue)

K2-663
100 pcs. with indicator.
made of plastic



SECURITY SEAL (blue)

K2-664
100 pcs. made of plastic



SECURITY SEAL (blue)

K2-661
100 pcs. with number.
made of plastic



SECURITY SEAL (red)

K2-662
100 pcs. with indicator.
made of plastic



SECURITY SEAL (red)

K2-660
100 pcs. made of plastic



SECURITY SEAL (red)

K2-665
100 pcs. with number.
made of plastic



DISPOSAL STAMP

K2-717
40 x 35 x 9 mm
Requested number or letter must be informed during the order.



LABEL HOLDER


K2-716
50 x 35 x 13 mm
Holder for identification labels

(-)




The Barrier Model is a mechanical filtering systems. The valves of Barrier Model Systems have three different positions during sterilization cycles such as relax position, before sterilization and after sterilization position. It can be used at least 5000 sterilization cycle.


1 Before sterilization cycle the valves of the Barrier Model Container are closed and in relax position.




2 At the beginning of the sterilization cycle, by the effect of vacuum the valves are opened completely and dry air goes out of the container.



3 During the sterilization cycle the valves are opened with effect of pressure in the reverse way totally opposite of the way moved by the effect of vacuum.



4 At the end of the sterilization cycle the autoclave (sterilizer) equates the pressure level and the valves opened the way by effect of vacuum.



5 After sterilization cycle the valves of the Barrier Model Container are closed and in relax position.

