

## Ophthalmological operating table – 2075 AR-EL



**Model 2075** has an electrically controlled vertical column height and the angle of inclination.

The base is equipped with wheels lock. The table features a seamless, antibacterial upholstery made of artificial leather. Metal parts are covered with ESD paint.

### Technical Details:

Table Height: 520mm - 880 mm  
Table Length: 2200mm  
Table Width: 550mm  
Table Base: 500mm x 980mm  
Head Section: 215mm x 235mm  
Shoulder Section: 480mm x 330mm  
Bed Section: 550mm x 1500 mm  
Lifting Capacity: 150 kg  
Motor: 2 Pieces  
Command Module: Foot switch

### Functional Details:

Head Section: Upwards: +20°  
Flexible polyurethane headpads bends all directions  
Shoulder Section: 0° ~ +35°  
Trendelenburg: 0° ~ +25°

### Standard Accessories:

Headrest frame (stainless steel), Anaesthetize frame (stainless steel), Rails&Armrests(2pcs)  
Optional Accessories: Locable Casters

### Packing Details:

Net weight: 157 Kg  
Gross weight: 195 Kg  
Packing Dimensions: 209×73×80cm

### Manufacturer:

Ar-El Berber ve Tibbi Ürünler San. ve Tic. Ltd. Şti.

