

2076-1 PROCTOLOGY TABLE

with 2 Motors Electrically adjustable height and trendelenburg – gas spring adjustable legrest



Technical Details:

Table Height : 690mm - 1050mm

Table Length : 1880mm

Table Width : 550mm

Table Base : 500mm x 980mm

Bed Section : 1000mm x 550mm

Leg Section : 500mm x 550mm

Knee Section : 300mm x 550mm

Lifting Capacity : 150 kg

Motor : 2 pieces

Command Module: Foot switch

Functional Details:

Trendelenburg : 0° ~ +23°

Leg Section : Downwards form the horizontal : -90°

Knee Section : Adjustable and movable up and down between 0 and 200mm.

Stomach support apparatus : Can easily assemble with both side rails and movable up and down.

Upholstery : Anti Bacterial , Artificial Leather

Standard Accessories:

Stomach support apparatus (1pc)

Optional Accessories : Lockable Casters

Packing Details:

Net weight : 160 Kg

Gross weight : 185 Kg

Packing Dimensions : 152×77×99cm



Manufacturer:

Ar-El Berber ve Tibbi Ürünler San. ve Tic. Ltd. Şti.

Distributor: **VIDEOMED ZAKŁAD ELEKTRONICZNY**

ul. Klonowa 18 58-310 Szczawno-Zdrój

tel. +48 74 843 81 09 fax +48 74 840 17 33

web: <http://www.videomed.eu> e-mail: handlowy@videomed.eu