

## Operating table – 2079-1 AR-EL

electrically adjustable height



### Parameters:

Adjustable height: 650mm - 1020mm  
Length: 1920mm  
Width: 550mm  
Base: 500mm x 980mm  
Headrest: 550mm x 260mm  
Back: 550mm x 500mm  
Seat: 550mm x 455mm  
Manhole: 550mm x 650mm  
Load capacity: 150 kg  
Engines: 1 piece  
Control: Foot

### Adjustments:

Headrest angle of inclination in the range of 0° ~ 80°	Removable, adjustable height 0-80 mm, adjustable angle of -90° ~ +10°
Backrest	angle of inclination in the range of 0° ~ 80°
Trendelenburg / Reverse Trendelenburg	22°/22°
Footrest	downward declination angle - 90°

### Standard accessories:

Drip hanger  
Anesthetic fraame  
Rails on both sides of the table  
Kneerests (2 pcs)

### Optional accessories:

Lockable wheels



### Producent:

Ar-El Berber ve Tibbi Ürünler San. ve Tic. Ltd. Şti.