

## Electrically adjustable gynecologic chair – 2082-1

with 3 Motors Electrically adjustable height – backrest – seat and trendelenburg



### Technical Details:

Table Height : 640mm - 990 mm  
Table Length : 1920mm (with foot extension)  
Table Width : 550mm  
Table Base : 500mm x 1000mm  
Seat Section : 550mm x 400 mm  
Back Section : 550mm x 600mm  
Head Section : 550mm x 280mm  
Foot extension : 550mm x 600mm  
Lifting Capacity : 150 kg  
Motor : 3 pieces  
Command Module : Foot switch

### Functional Details:

Head Section : Detachable , Adjustable height between 0 - 80mm  
Back Section : 0° ~ 70°  
Trendelenburg / Reverse trendelenburg : -5° ~ +20°  
Foot extension : Detachable

### Standard Accessories :

Foot Extension (1pc), Knee crutches: (2pcs), Handle (2pcs), Waste tray (1pc), Rails&Armrest (2pcs), Table covering (1pc), Movable head cushion (1pc), Paper roll holder (1pc), Lockable Casters (3pcs)

### Packing Details:

Net weight : 190 Kg  
Gross weight : 215 Kg  
Packing Dimensions : 152×77×99cm

### Manufacturer:

Ar-El Berber ve Tıbbi Ürünler San. ve Tic. Ltd. Şti.