

Electrically adjustable gynecologic chair – 2088

electrically adjustable height



Technical Details:

Seat Height : 620mm - 920 mm
Table Length : 1950mm (with extension)
Table Width : 640mm
Seat Section : 550mm x 360mm
Back Section : 580mm x 965mm
Foot extension : 550mm x 600mm
Lifting Capacity : 200 Kg
Motor : 3 pieces
Command Module : Foot switch

Functional Details:

Back Section : 0° ~ +70°
Trendelenburg / Reverse trendelenburg: -5° ~ +25 °
Foot extension : Detachable

Standard Accessories :

Foot extension (1 pc), Knee crutches (2 pcs),
Handle (2 pcs), Waste tray (1 pc), Table covering(1pc)
Paper roll holder (1 pc), Accessory rail (2 pcs)

Packing Details:

Net weight: 128 Kg
Gross weight :153 Kg
Packing Dimensions : 152×77×99cm



Manufacturer:

Ar-El Berber ve Tıbbi Ürünler San. ve Tic. Ltd. Şti.