

Examline 850 Blood Donor Chair

Ref: M2Z 000 0 0



Bıçakcılar Examline 850 Blood Donor Chair is designed and manufactured to provide years of trouble-free service in a variety of applications such as: Blood Banks, Dialysis, Pheresis and Plasma Centers, Hospitals and Clinics.

Ergonomically designed and constructed from the finest materials.

Adjustable armrests provide both patient and operator comfort. A linear actuator motor provides smooth positioning. The convenient hand pendant, which is easily attached to either side of the head of the chair, controls the motor allowing easy positioning from "head high-foot low" to "head low-foot high" or any suitable position in between.

Technical Details

- Seat height:** 500 mm
- Head height:** 1050 - 860 mm
- Foot height:** 650 - 840 mm
- Chair length:** 1470 mm
- Chair width:** 650 mm
- Chair width (w/armrest sections):** 980 mm
- Voltage:** 230V, 50/60Hz

Packing Details

- Net weight:** 52 kg
- Gross weight:** 133 kg
- Packing dimensions:** 900 x 1700 x 1160 mm

Positions

Head and foot heights:
Adjustable



Standard Accessories

- Armrests (2 ea)
- Legrest cover (1 ea)



Armrests - Adjustable, detachable

Ref: M2Z 065 4 0 - Right
Ref: M2Z 064 4 0 - Left



Legrest cover

Ref: M2Z 078 5 0

Accessories

- I.V. stand



I.V. stand

Ref: M2Z 100 5 0