



Colposcopes AC-1 series

Available models:

- AC-1100 halogen light-source (6V / 20W), one magnification of 9x
- AC-1100 LED light-source (5,9V / 5W), one magnification 9x
- AC-1300 halogen light-source (6V / 20W), three magnifications levels of 3.7x, 7x, 15x
- AC-1310 LED light-source (5,9V / 5W), three magnifications levels: 3.7x, 7x, 15x
- AC-1320 LED light-source (5,9V / 5W), three magnifications levels: 3.7x, 7x, 15x; binocular on 45° angle

Thanks to it's lightweight and mounted wheels, colposcopes AC-1 series can be easily moved.

Specifications:

- magnification: AC-1110: 9x AC-1310/1320: 3.75x, 7x, 15x focal AC-1110: 316mm AC-1310/1320: 300mm field of view: AC-1110: 22mm AC-1310/1320:30.93mm, 15.50mm, 7.74mm depth of field: AC-1110: > 4mm AC-1310: 18mm, 4mm, 2mm AC-1320: 22mm, 4mm, 3mm light intensity • AC-1110: > 15000 lux AC-1310: > 28000 lux AC-1320: > 10000 lux individually adjustable glasses 12,5 • two measuring circles, built in to the right eyepiece
- continuously adjustable brightness
- green light filter
- easily replaceable bulb
- height adjustment: 838mm ~ 1194mm
- Universal Teflon® ball joint that does not require pressure control and allows for smooth rotation of 360°
- the provided product is ready to use out of the box (no need to call a technician for installation)





Dwa okręgi pomiarowe w modelach 1310 i 1320

Manufacturer:

Alltion (Wuzhou) Co., Ltd. Suite B, 11/F, Lijing Tower, Cathay Plaza, NO.18, Zhongshan Road, Wuzhou, Guangxi, China

Distributor: **VIDEOMED ZAKŁAD ELEKTRONICZNY** ul. Klonowa 18 58-310 Szczawno-Zdrój tel. +48 74 843 81 09 fax +48 74 840 17 33 web: http://www.videomed.eu e-mail: handlowy@videomed.eu

