

Camera C3.multi



LUT

Camera C3.multi

LED ON/OFF indicator, white balance and window function. Preset function with three selectable basic setups, modifiable by OSD-steering, brightness adjustment with OSD indicator. Camera head with two push buttons: remote, gain.

Camera controller

Outputs: 1 x RGB Sub-D, 2 x S-Video (Y/C), 1 x Video (CVBS),
1 x 3.5 mm cinch for remote control

Camera head

Image sensor: 3 x 1/3" CCD

Aspect ratio: 4:3

Cable length: 3 m

Endocoupler: f = 23 mm

System

Image resolution: 750x430 (PAL), 750x350 (NTSC) lines

Medical device directive 93/42/EEC Class 1, EN 60601-1; BF/CE mark

100 - 250 V AC, 50/60 Hz, 25 VA

Dimensions / Weight

Camera controller: (W/H/D) 360x88x320 mm

Camera head: (L/Ø) 165x50 mm

Weight camera controller: approx. 4.0 kg

Order number

05.9090s.m1

The logo for LUT GmbH, featuring the letters 'LUT' in a large, bold, white, stylized font with a thick outline, set against a dark green background.

LUT GmbH

Robert-Bosch-Strasse 3

D - 79211 Denzlingen

www.lut-endoscopy.com

Tel.: +49 (0) 76 66/9 00 69-0

Fax: +49 (0) 76 66/9 00 69-480

info@lut-gmbh.de

German Endoscopy Manufacturer