

CHEESE™



Mini Dental Diode Laser Systems

4W/7W/10W

Dentist's Good Assistant!

+ Features

- User-friendly interface
- Customized programs
- 810nm/980nm for choice
- Battery operated



Wuhan Gigaa Optronics Technology Co.,Ltd.

www.gigaalaser.com

Mini Dental Diode Laser Systems

Dentist's Good Assistant



CHEESE™ dental laser systems, the top class soft tissue laser in the world provided by GIGAA laser. Dentists like our dental laser and make the patient smile, sound like "CHEESE!"

Mini dental laser CHEESE™ could perform a wide range of soft-tissue surgical, periodontal, and endodontic procedures. There are 17 presets treatment programs in the touch screen software and two wavelengths of 810nm/980nm for your choice, the innovative battery-operated design will provide out door visiting much convenience.

CHEESE™ laser is an ideal choice for your dental clinic and will make it attractive.

Clinical Applications

+ Laser Periodontal

- Preprocedural decontamination
- Laser decontamination
- Coagulation
- Gingivectomy/curettage

+ Laser Aesthetics & Surgery

- Abscess
- Aphthae/ulcers
- Fibroma
- Frenectomy
- Gingivectomy
- Hemostasis/coagulation
- Implant recovery
- Leukoplakia
- Crown lengthening
- Endodontics

+ Laser Endodontics

- Root canal sterilization
- Desensitization

+ Laser Teeth Whitening(7W only)



Surger Handpiece With Disposable Tips



Bleaching Handpiece



Laser Type	GaAlAs diode Laser		
Model	DEN4A/B	DEN7A/B	DEN10B/C
Wavelength	810nm/980nm	810nm/980nm	980nm/940nm
Maximum Power	4W	7W	10W
Operation Mode	CW, Single or Repeat Pulse		
Pulse Duration	1ms-1s		
Repetition Rate	0.5Hz-0.5KHz		
Fibers	Fibers Of 200µm ,400µm And 600µm With SMA905 Connector		
Pilot Beam	Red Diode Laser of 635nm, Power<3mW		
Control Mode	Micro Touch Key and Touch Screen		
Voltage/Battery Power	Adapter:AC110/240V~47/63 Hz,1.25-0.5A/Battery:DC 25.7/3A		
Dimensions	120(W)*200(L)*155(H) mm		
Weight	1.8Kg		

Clinic cases

Please visit our website www.gigaalaser.com for more videos.

Accessories




Dental Handpiece



Fiber



Safety Goggle



Fiber Cutter



Foot Switch



Fiber Holder



Battery



ADD: B8-A5, Hi-Tech Medical Device Industrial Park, #818 Gaoxin Avenue,
Wuhan 430206, China
Ph: 86-27-67848871 67848872 Fax: 86-27-67848873
<http://www.gigaalaser.com> Email: info@gigaalaser.com

