

SURGICAL SUCTION UNIT Quality Speaks Itself DF-650

- 1. The outer shall of the instrument adopted by ABS material. Clear simple lines and nice looking appearance create a new enjoyment for users. Humanized design with outstanding technology impressed users in big hospitals, medium/small clinics, nursing homes and home cares.
- 2. Adhering to the concept of environmental protection, ecological protection and econmical benefits makes each item of material be recycles as possible as it can to reduce waste and protect natural resource.

Functions:

- Digital Vacuum Monitor is used to show the suction power, so users can easily read it.
- The suction tube directly connects to the suction bottle, so users can easily wash it. It's sanitary & safe. 3000 c.c. bottle
- The handle attached to the suction bottle makes users easily lift and wash it.
- The foot switch adopts air pressure type, so it's anti-moisture, insulating & safe.
- 360° Individual brake wheel features easy movement and high safety.
- The fastening rack installed at the side makes the suction tube be easily put.

Specifications

Vacuum Range

250 W 1/3 HP

10~720 mm/hg

3000c.c. x 2 pcs

Air pressure type

L45xW35xH84 cm

80 liter/min

Max. HP

Max. Flow

Foot control

Machine Size

Machine Weight 19 KGS

Capacity

The Instrument plane:

- Digital vacuum monitor
- Regulator valve
- Surgical tube holder
- Surgical suction adapter
- 6 Power switch (On/Off)
- Filter (L)

The Side:

Power cord

3 Circuit breaker

Adapter of foot switch

4 All casters with locking breaks 5 Foot switch (Air pressure type)

- 8 Floating control (L)
- Floating control (S)













OIL-LESS VACUUM PUMP (90-026) Twin cylinder piston pump design.



Standard accessories

1. Suction tube (S) 180cm 1 pc

2. Suction catheter 1 pc

3. Foot switch 1 pc

4. Filter2 pcs

Humanizing design: Low noise, Wear resistance.

Compact structure: Low power consumption, long life.

Fully oil-less: Non-pollution, free maintenance

Voltage	120V	240V
Current	2.4A	1.2A
Capacitor	26uf	13uf
R.P.M.	1720	
Poles	4	
Output	250W 1/3HP	
Max. flow rate	80 liter/min	
Max. suction	700~720 mm/hg	
Weight	9.0 kgs	
Thermal protector, auto reset		

Pump specification



1 Fix cover 2 Depositing box 3 Fan

