

Transportable suction unit DF-650B



Suction unit DF-650B has been designed specifically for plastic surgery, especially for liposuction. DF-650B most characteristic feature is its very high flow, up to 112 l/min.

The unit comes with two reusable bottles, each with a volume of 3000ml. Bottles are designed for sterilization. After filling the first bottle the device automatically switches to the other one.

DF-650B has safeguards to protect the vacuum pump from being flooded. An additional protection is a small tank with an anti-overflow valve.

The device is very easy to move thanks to the 4 wheels with a turn radius of 360°. Each of the wheels has a brake.

The suction unit has a large gauge which allows the user to easily read the values. ABS plastic is used for the housing of DF-650B making it very easy to clean and gives it a slightly shiny appearance.

Specifications:

- Power: 250W
- Maximum flow: > 112 l / min
- Maximum pressure: >720 mmHg (95kPa)
- Continues negative pressure control: 0-720 mmHg
- Wheels: 4 with breaks
- Dimensions: 450×400×960
- Weight: 32 kg

Standard accessories:

- 2x Suction tube (S) 180 cm
- 1x Suction tube (L) 180 cm
- 1x Suction Catheter
- 1x Footswitch
- 2x Filter (L)
- Bottles: 2×3 l (2×2 l – optional)

Manufacturer:

Doctor's Friend Medical Instrument Co., Ltd.
No.28, Lane 13, Sec. 3, Zhong Shan Rd., Tanzi Shiang, Taichung County, 42756, Taiwan (R.O.C.)