



Portable suction units – DF-760 & DF-760A

DF-760



FEATURES:

- Whole A.B.S plastic injection made outer shell. Design with smart appearance, light weight and easy to carry. Low noise and powerful suction function which provides a better quality care for users.
- Regulator was designed to adjust the pressure which could help suction function becomes more efficient. The operation plate was made by PC which can prevent from scratch and discoloration.
- While the bottle is full, the level float control will stop the suction function automatically which could prevent the pollution against the pump.
- DF-760 is the necessary item of hospital,

DF-760A



clinics, pediatrics and home care.

STANDARD ACCESSORIES:

- 1x Suction tube 180cm
- 1x Suction catheter
- 2x Filter(S)

SPECIFICATIONS:

	DF-760	DF-760A
Horse power	88W	88W
Flow rate	35 liter/min	35 liter/min
Vacuum pressure	65~70 cm/Hg	65~70 cm/Hg
Capacity	1000c.c.	1000c.c. x 2pcs
Safe device	Level float control	Level float control
N.W. / G.W.:	6.3 / 7.0 kg	6.9 / 8.0 kg

Doctor's Friend Medical Instrument Co., Ltd. No.28, Lane 13, Sec. 3, Zhong Shan Rd., Tanzi Shiang, Taichung County, 42756, Taiwan (R.O.C.)