

## 2077-4 DIALYSIS CHAIR

with 4 Motors Electrically adjustable height – backrest – trendelenburg – legrest



### Technical Details:

Seat height : 520mm - 950mm  
Chair Length : 2130mm  
Chair Width : 700mm, (with Armrests: 1060mm)  
Headrest : 700mm x 300mm ( adjustable 100mm )  
Seat Section : 700mm x 540mm  
Back Section : 700mm x 600mm  
Leg Section : 700mm x 620mm  
Foot support : 200mm x 650mm ( adjustable 250mm )  
Lifting Capacity : 200 kg  
Command Module : Hand switch

Motor : 4 pieces

### Functional Details:

Back Section :  $0^{\circ} \sim +70^{\circ}$   
Leg Section :  $0^{\circ} \sim +30^{\circ}$   
Trendelenburg / Reverse trendelenburg :  $-5^{\circ} \sim +15^{\circ}$   
Headrest: Adjustable upwards- downwards 100mm  
Armrest : Adjustable upwards- downwards.  
Serum hanger: Can be use both left and right side  
Foot support: Adjustable upwards- downwards 250mm

### Standart Accessories:

Serum hanger, Movable head cushion.

### Packing Details:

Net weight : 115 Kg  
Gross weight: 145 Kg  
Packing Dimensions : 185x110x70cm

Warranty : 2 year



## 2077 DIALYSIS & BLOOD DONOR CHAIR

with 1 motor electrically adjustable backrest or Gas spring adjustable backrest – seat and Legrest moves together with backrest



### Technical Details:

Seat height	: 550 mm
Chair Length	: 1730mm
Chair Width	: 800mm ( including tripod: 1300mm)
Seat Section	: 530mm x 500mm
Back Section	: 530mm x 810mm
Leg Section	: 530mm x 420mm
Lifting Capacity	: 150 kg
Motor	: 1 piece or manual (optional)
Command Module	: Hand switch or manual (optional)

### Functional Details:

Head Section	: Separable
Back Section	: $-5^{\circ} \sim +80^{\circ}$
Seat Section	: Moves together with backrest ( Synchronized )
Armrest	: Fixed
Dialysis tripods	: Adjustable parallel or vertical to armrest
Leg Section	: Moves together with backrest ( Synchronized ) or independent motion.

### Standart Accessories:

Serum hanger, Examination tissue paper , Movable head cushion

### Packing Details:

Net weight	: 70 Kg
Gross weight	: 90 Kg
Packing Dimensions	: 102x82x155cm

**Warranty : 2 years**



## 2008 DIALYSIS & BLOOD DONOR CHAIR with 2 Motors Electrically adjustable Backrest and Legrest – fixed Seat



### Technical Details:

Seat height	: 500 mm
Chair Length	: 1600mm
Chair Width	: 600mm
Chair Width	: 800mm (with Armrests)
Seat Section	: 600mm x 450mm
Back Section	: 600mm x 770mm
Leg Section	: 560mm x 360mm
Lifting Capacity	: 150 kg
Command Module	: Hand switch

### Functional Details:

Back Section	: -8° ~ +85°
Seat Section	: Fixed
Armrest	: Adjustable upwards- downwards & left- right & forward - back and Swivel
Leg Section	: Foldable, 0° ~ 90°

### Standart Accessories:

Serum hanger, Movable head cushion

### Packing Details:

Net weight	: 48 Kg
Gross weight	: 68 Kg
Packing Dimensions	: 135x102x82cm

**Warranty : 2 years**



## 2007 BLOOD DONOR CHAIR

Gas spring adjustable backrest – Legrest moves together with backrest – fixed seat



### Technical Details:

Seat height	: 430 mm
Chair Length	: 1620mm
Chair Width	: 760mm
Seat Section	: 550mm x 410mm
Back Section	: 550mm x 760mm
Leg Section	: 500mm x 300mm
Lifting Capacity	: 150 kg
Command Module	: Manual

### Functional Details:

Back Section	: -5° ~ +75°
Seat Section	: Fixed
Armrest	: Adjustable height and Swivel
Leg Section	: Moves together with backrest ( Synchronized )

### Standart Accessories:

Serum hanger (1pc), Armrests (2pcs)

### Packing Details:

Net weight	: 48 Kg
Gross weight	: 68 Kg
Packing Dimensions	: 135x102x82cm

**Warranty : 2 years**

## 2007-1 VISITOR CHAIR

Gas spring adjustable backrest – Legrest moves together with backrest – fixed seat



### Technical Details:

Seat height	: 470 mm
Chair Length	: 1700mm
Chair Width	: 680mm
Seat Section	: 580mm x 470mm
Back Section	: 580mm x 900mm
Leg Section	: 550mm x 300mm
Lifting Capacity	: 180 kg
Command Module	: Manual

### Functional Details:

Back Section	: 0° ~ +85°
Seat Section	: Fixed
Armrest	: Fixed
Leg Section	: Moves together with backrest between 0° to 90°
Upholstery	: Anti Bacterial , Artificial Leather

### Packing Details:

Net weight	: 48 Kg
Gross weight	: 68 Kg
Packing Dimensions	: 172x70x 70cm

**Warranty : 2 years**