

1001 DOCTOR STOOL



1002 DOCTOR STOOL



1001 DOCTOR STOOL

Stool height : minimum : 520mm , maximum 670mm

Stool backrest : up to 9 cm movements.(up and down)

Stool backrest : Adjustable beetween 90° to 115°

Backrest angle:Adjustable backrest angle with circular buttons, which located under the seat section.

Right side button adjusts backrest height, Left side button adjusts backrest angle.

Seat height : Adjustable seat height with pushing circular metal foot treadle.

Pushing metal foot treadle while sitting on stool, the stool height gets lower.

Pushing metal foot treadle while stool is empty, the stool height gets upper.

Upholstery : Anti Bacterial , Artificial Leather

Packing Details:

Net weight :14 Kg

Gross weight :17 Kg

Packing Dimensions : 62x57x58cm

1002 DOCTOR STOOL

Stool height : minimum : 500mm , maximum 650mm

Stool cushion : diameter : 38 cm

Buttom Circular Metal : diameter : 60 cm

Seat height : Adjustable seat height with pushing circular metal foot treadle.

Pushing metal foot treadle while sitting on stool, the stool height gets lower.

Pushing metal foot treadle while stool is empty, the stool height gets upper.

Upholstery : Anti Bacterial , Artificial Leather

Packing Details:

Net weight : 11 Kg

Gross weight :13 Kg

Packing Dimensions : 62x57x58cm

1003 DOCTOR STOOL



1004 DOCTOR STOOL



1003 DOCTOR STOOL (Hydraulic)

Stool height : minimum : 510mm , maximum 670mm

Stool backrest : up to 9 cm movements.(up and down)

Stool backrest : Adjustable between 90° to 115°

Backrest angle:Adjustable backrest angle with circular buttons, which located under the seat section.

Right side button adjusts backrest height, Left side button adjusts backrest angle.

Seat height : Adjustable seat height with pushing circular metal foot treadle.

Pushing metal foot treadle while sitting on stool, the stool height gets lower.

Pushing metal foot treadle while stool is empty, the stool height gets upper.

Upholstery : Anti Bacterial , Artificial Leather

Packing Details:

Net weight :14 Kg

Gross weight :17 Kg

Packing Dimensions : 62x57x58cm

1004 DOCTOR STOOL(Hydraulic)

Stool height : minimum : 450mm , maximum 610mm

Stool cushion : diameter : 38 cm

Buttom Circular Metal : diameter : 60 cm

Seat height : Adjustable seat height with pushing circular metal foot treadle.

Pushing metal foot treadle while sitting on stool, the stool height gets lower.

Pushing metal foot treadle while stool is empty, the stool height gets upper.

Upholstery : Anti Bacterial , Artificial Leather

Packing Details:

Net weight : 11 Kg

Gross weight :13 Kg

Packing Dimensions : 62x57x58cm