

2087 GYNECOLOGICAL OBSTETRIC DELIVERY TABLE

with 4 Motors Electrically adjustable height – backrest – trendelenburg – Knee Crutches

with 3 Motors Electrically adjustable height – backrest – trendelenburg manual Knee Crutches



Technical Details:

Seat Height : 640mm - 940 mm

Table Length : 2010mm

Table Width : 640mm

Seat Section : 640mm x 360mm

Back Section : 640mm x 1000mm

Footrest Section : 640mm x 600mm

Curved Leg Section : Min space:850mm

Lifting Capacity : 200 Kg

Motor : 4 or 3 pieces

Command Module : Foot switch

Functional Details:

Back Section : 0° ~ +50°

Trendelenburg : -5° ~ +25°

Foot Extension : Detachable

Curved Legrest : Adjustable upwards-downwards & left- right and Detachable

Standart Accessories :

Foot Extension (1 pc), Curved Legrest (2 pcs),

Knee crutches (2 pcs), Handle (2 pcs),

Waste tray (1 pc), Paper roll holder (1 pc),

Accessory hanger (1 pc).

Optional Accessories :

Foot Step (2 pcs)

Packing Details:

Net weight: 139 Kg

Gross weight :169 Kg

Packing Dimensions : 152 x 77 x 99 cm

Packing Type : Full protect wooden case.

Warranty : 2 year

