



Portable LED light source LL020



Ultra bright LED light source. High light efficiency thenks to an optimized endoscope connection. Control of light intensity, secure grip thanks to the grooved surface. The light source is small and handy. Compatible with all endoscopes from the LUT company. Comes with battery and a charger.

Lamp:

Type: high power LED Illumination strength: 20 000 lx Color temperature: 5500°K Work time: ~45 min (at the highers brightness) after every charge up

Control:

Light intensity control Optic fiber socket for Storz adapter Battery: Li-Ion, 950 mAh, rechargeable Charging time: ~2 h

Dimensions / Weight: Dimensions / Ø: 89 mm / 27 mm Weight: ~90g (with battery)

Light source LED LL020, z battery and charger Additional battery Li-Ion, 950 mAh 05.0755led 05.0755led.bat

Distributor: VIDEOMED ZAKŁAD ELEKTRONICZNY ul. Klonowa 18 58-310 Szczawno-Zdrój tel. +48 74 843 81 09 fax +48 74 840 17 33 web: http://www.videomed.eu e-mail: handlowy@videomed.eu