

2042-1 ENT CHAIR

Hydraulic adjustable height – gas spring adjustable backrest



Technical Details:

Chair Height : 560mm - 700 mm

Chair Width : 630mm

Chair Base : 610

Seat Section : 500mm x 420 mm

Back Section : 500mm x 650mm

Head Section : 210mm x 210mm

Lifting Capacity : 150 kg

Command Module : Hydraulic foot pedal

Functional Details:

Head Section : Detachable , Adjustable front and back, adjustable height between 0 -150mm.

Back Section : Adjustable between 50° ~ 75°

Seat Section : Swivel 360°

The act of rotation could be fixed with pushing hydraulic foot pedal upwards.

Packing Details:

Net weight : 55 Kg

Gross weight : 71 Kg

Packing Dimensions : 70x80x126cm

Warranty : 2 years