

2060 ENT CHAIR

with 2 Motors Electrically adjustable height and backrest



Technical Details:

Chair Height : 480mm - 880 mm
Chair Length : 1950mm (without headrest: 1650mm)
Chair Width : 620mm (with armrests: 700mm)
Chair Base : 500mm x 700mm
Seat Section : 500mm x 580 mm
Back Section : 500mm x 540mm
Head Section : 210mm x 210mm
Foot extension section : 540mm x 610mm
Foot Step : 350mm x 230mm
Lifting Capacity : 200 kg
Motor : 2 pieces
Command Module : Hand switch

Functional Details:

Head Section : Detachable , Adjustable front and back, adjustable height between 0 -150mm.
Back Section : 0° ~ +90°
Foot extension : Detachable
Armrest : Detachable

Packing Details:

Net weight : 97 Kg
Gross weight : 122 Kg
Packing Dimensions : 152x77x99cm

Warranty : 2 years