

## 2061-2 ENT CHAIR ( SWIVEL ) with 2 Motors Electrically adjustable height and backrest Mechanical adjustable swivel



### Technical Details:

Seat Height : 590mm - 990 mm  
Chair Length : 1900mm, ( without headrest: 1600mm )  
Chair Width : 650mm, ( without armrests: 550mm )  
Chair Base : 500mm x 700mm  
Seat Section : 550mm x 450 mm  
Back Section : 500mm x 540mm  
Head Section : 210mm x 210mm  
Leg Section : 510mm x 260 mm  
Foot Step : 435mm x 250mm  
Motor : 2 pieces  
Lifting Capacity : 200 kg  
Command Module : Hand switch

### Functional Details:

Rotation : 180° left and 180° right and locable.  
Head Section : Detachable , Adjustable front and back, adjustable height between 0 -150mm.  
Back Section : 0° ~ +85°  
Foot step : Foldable

### Packing Details:

Net weight : 110 Kg  
Gross weight : 135 Kg  
Packing Dimensions : 152x77x99cm

**Warranty : 2 years**