

2101 ENT CHAIR

with 1 Motor Electrically adjustable height – gas spring adjustable backrest



Technical Details:

Chair Height : 560mm - 930 mm
Chair Length : 1550mm
Chair Width : 600mm
Chair Base : 500mm x 980mm
Seat Section : 550mm x 490 mm
Back Section : 500mm x 540mm
Head Section : 210mm x 210mm
Foot step Section : 320mm x 220mm
Lifting Capacity : 150 kg
Motor : 1 piece
Command Module : Hand switch

Functional Details:

Head Section : Detachable , Adjustable front and back, adjustable height between 0 -150mm.
Back Section : Adjustable between 40° ~ 85°
Armrest : Adjustable vertical 90°
Foot step : Detachable

Packing Details:

Net weight : 135 Kg
Gross weight : 160 Kg
Packing Dimensions : 152x77x99cm

Warranty : 2 years