

2061-2 ENT CHAIR (SWIVEL) with 2 Motors Electrically adjustable height and backrest Mechanical adjustable swivel



Technical Details:

Seat Height : 590mm - 990 mm
Chair Length : 1900mm, (without headrest: 1600mm)
Chair Width : 650mm, (without armrests: 550mm)
Chair Base : 500mm x 700mm
Seat Section : 550mm x 450 mm
Back Section : 500mm x 540mm
Head Section : 210mm x 210mm
Leg Section : 510mm x 260 mm
Foot Step : 435mm x 250mm
Motor : 2 pieces
Lifting Capacity : 200 kg
Command Module : Hand switch

Functional Details:

Rotation : 180° left and 180° right and locable.
Head Section : Detachable , Adjustable front and back, adjustable height between 0 -150mm.
Back Section : 0° ~ +85°
Foot step : Foldable

Packing Details:

Net weight : 110 Kg
Gross weight : 135 Kg
Packing Dimensions : 152x77x99cm

Warranty : 2 years

2060 ENT CHAIR

with 2 Motors Electrically adjustable height and backrest



Technical Details:

Chair Height : 480mm - 880 mm
Chair Length : 1950mm (without headrest: 1650mm)
Chair Width : 620mm (with armrests: 700mm)
Chair Base : 500mm x 700mm
Seat Section : 500mm x 580 mm
Back Section : 500mm x 540mm
Head Section : 210mm x 210mm
Foot extension section : 540mm x 610mm
Foot Step : 350mm x 230mm
Lifting Capacity : 200 kg
Motor : 2 pieces
Command Module : Hand switch

Functional Details:

Head Section : Detachable , Adjustable front and back, adjustable height between 0 -150mm.
Back Section : 0° ~ +90°
Foot extension : Detachable
Armrest : Detachable

Packing Details:

Net weight : 97 Kg
Gross weight : 122 Kg
Packing Dimensions : 152x77x99cm

Warranty : 2 years

2101 ENT CHAIR

with 1 Motor Electrically adjustable height – gas spring adjustable backrest



Technical Details:

Chair Height : 560mm - 930 mm
Chair Length : 1550mm
Chair Width : 600mm
Chair Base : 500mm x 980mm
Seat Section : 550mm x 490 mm
Back Section : 500mm x 540mm
Head Section : 210mm x 210mm
Foot step Section : 320mm x 220mm
Lifting Capacity : 150 kg
Motor : 1 piece
Command Module : Hand switch

Functional Details:

Head Section : Detachable , Adjustable front and back, adjustable height between 0 -150mm.
Back Section : Adjustable between 40° ~ 85°
Armrest : Adjustable vertical 90°
Foot step : Detachable

Packing Details:

Net weight : 135 Kg
Gross weight : 160 Kg
Packing Dimensions : 152x77x99cm

Warranty : 2 years

2042-1 ENT CHAIR

Hydraulic adjustable height – gas spring adjustable backrest



Technical Details:

Chair Height : 560mm - 700 mm
Chair Width : 630mm
Chair Base : 610
Seat Section : 500mm x 420 mm
Back Section : 500mm x 650mm
Head Section : 210mm x 210mm
Lifting Capacity : 150 kg
Command Module : Hydraulic foot pedal

Functional Details:

Head Section : Detachable , Adjustable front and back, adjustable height between 0 -150mm.
Back Section : Adjustable between 50° ~ 75°
Seat Section : Swivel 360°
The act of rotation could be fixed with pushing hydraulic foot pedal upwards.

Packing Details:

Net weight : 55 Kg
Gross weight : 71 Kg
Packing Dimensions : 70x80x126cm

Warranty : 2 years