

Air caloriser **KALORistar Arctic** with active cooling



KALORistar Arctic is a powerful device for thermal stimulation of the middle ear. It has an active cooling system, so it is possible to stimulate with cold air even on hot days. Recommended for private offices and clinics. Arctic can be controlled from a PC computer using eVNG software by BioMed Jena GmbH.

Specification:

- Power supply: 230V + / - 5%, 50Hz
- Power consumption: max. 250 mA
- Power consumption: max. 60W
- Stimulus duration: 10s ~ 210s
- Air temperature: 20° C ~ 47° C (max. 8° C below ambient temperature)
- Air flow: 5.0 l / min. + / - 10%
- Dimensions: 115 x 307 x 257 mm
- Weight: 3.5 kg

Safety:

- Safety class (EN 60601-1): I
- Degree of safety; Type B
- Type of protection: IPX0
- IIA Classification: according to Encl. IX Guidelines EG 93/42/EEC
- CE Mark: CE 0124
- Applied standards: EN 60601-1:1990 + A1: 1993 + A2: 1995, EN 60601-1-2:2001 (EMV/EMC)

Manufacturer:

BioMed Jena GmbH

Lutherstr. 148, 07743 Jena, Germany