

STERILE CONTAINER SYSTEMS A1 MODEL



1/1 FULL SIZE



3/4 QUARTER SIZE



1/2 HALF SIZE

Aygün A1 model sterilization container systems basically have similar features with standard model sterilization container systems. It has been designed by restricting some features of standard model which will not risk any of sterilization necessities. When compared to standard model sterilization container systems there is no colorful lid option in this model and coloring parts are carried out to label plates. It can be supplied with five main colors except for black. Moreover the handles of this model are not silicone coated and this model cannot be supplied with perforated boxes and more than 200 mm deep boxes. Similar to other models, Aygün A1 model sterilization container systems are manufactured as three different length, width and height. All models and sizes have specific baskets and accessories which can be viewed in following parts.

STERILE CONTAINER SYSTEMS



A1 Model



» *1/1 Full Size*

1/1 FULL SIZE



W: 285 mm





















L: 600 mm



Lid (Perforated)



Box (Non Perforated)

h:110 mm	 K6-010E	 K6-011E	 K6-012E	 K6-013E	 K6-014E
h:135 mm	 K6-016E	 K6-017E	 K6-018E	 K6-019E	 K6-020E
h:155 mm	 K6-022E	 K6-023E	 K6-024E	 K6-025E	 K6-026E
h:205 mm	 K6-028E	 K6-029E	 K6-030E	 K6-031E	 K6-032E



STANDARD MODEL
29

BIO-BARRIER MODEL
63

Polymer's PPSU Lids
97

WINDOW MODEL
119

B₃ Barrier MODEL
139

PLASMA MODEL
173

A1 MODEL
185

Dental, Mini, Special
205

BASKETS
223


Autoclave Mechanism Emulating Indicator
243

ACCESSORIES
251

BOX (non perforated)




outside dimension



outside dimension



inside dimension



outside dimension	inside dimension	weight
600 x 272 h:100 mm	543 x 258 h:94 mm	1900 g
600 x 272 h:116 mm	543 x 258 h:110 mm	2096 g
600 x 272 h:138 mm	543 x 258 h:132 mm	2205 g
600 x 272 h:190 mm	543 x 258 h:184 mm	2648 g


PERFORATED LID



Dimension: 602 x 285 mm
Weight : 1300 g

BASKET MADE FROM STAINLESS STEEL SHEET

BASKET (with stand)



	code	outside dimension
	K3-740	540 x 253 h:40 mm
	K3-741	540 x 253 h:60 mm
	K3-742	540 x 253 h:80 mm
	K3-743	540 x 253 h:110 mm

➤ The baskets which are without lid and less than 44 mm in height do not include the handles while the ones with lid include the handles.


STAND



code	height
K3-767	20 mm




STAND



code	height
K3-768	10 mm



STAND



code	height
K3-766	5 mm




BASKET MADE FROM STAINLESS STEEL WIRE MESH




code	outside dimension
K2-460	540 x 245 h:30 mm
K2-461	540 x 245 h:50 mm
K2-462	540 x 245 h:70 mm
K2-463	540 x 245 h:100 mm

SILICONE MAT



K2-526
520 x 230 mm
K2-528
475 x 230 mm







SILICONE MESH MAT



K2-531
485 x 250 mm

SILICONE HANDLE




K2-710 
K2-711 
K2-712 
K2-713 
K2-714 
K2-715 

PAPER FILTER




K2-740
Single use 100 pcs
with indicator. Aygun mark.
Diameter: 184 mm

ALUMINIUM FILTER RETAINER




K2-722
Diameter: 186 mm







FILTER DISPENSER




K2-730

IDENTIFICATION LABELS (aluminium)



K2-700 
K2-701 
K2-702 
K2-703 
K2-704 
K2-705 


LABELS (paper)



K2-673
100 pcs. With indicator.

For small size and
special containers

LABELS (paper)



K2-672
100 pcs. With indicator.

For full and
Half size containers

SECURITY SEAL (blue)



K2-663
100 pcs. with indicator.
made of plastic

SECURITY SEAL (blue)



K2-664
100 pcs. made of plastic

SECURITY SEAL (blue)



K2-661
100 pcs. with number.
made of plastic

SECURITY SEAL (red)



K2-662
100 pcs. with indicator.
made of plastic

SECURITY SEAL (red)



K2-660
100 pcs. made of plastic

SECURITY SEAL (red)



K2-665
100 pcs. with number.
made of plastic

STANDARD MODEL	29
BIO-BARRIER MODEL	63
Polymers PPSU Lids	97
WINDOW MODEL	119
B ₃ Barrier MODEL	139
PLASMA MODEL	173
A1 MODEL	185
Dental, Mini, Special	205
BASKETS	223
Autolock Mechanism Emulating Indicator	243
ACCESSORIES	251

A1 MODEL

1/1 full size

