

STERILE CONTAINER SYSTEMS A1 MODEL



1/1 FULL SIZE



3/4 QUARTER SIZE



1/2 HALF SIZE

Aygün A1 model sterilization container systems basically have similar features with standard model sterilization container systems. It has been designed by restricting some features of standard model which will not risk any of sterilization necessities. When compared to standard model sterilization container systems there is no colorful lid option in this model and coloring parts are carried out to label plates. It can be supplied with five main colors except for black. Moreover the handles of this model are not silicone coated and this model cannot be supplied with perforated boxes and more than 200 mm deep boxes. Similar to other models, Aygün A1 model sterilization container systems are manufactured as three different length, width and height. All models and sizes have specific baskets and accessories which can be viewed in following parts.

STERILE CONTAINER SYSTEMS




A1 Model

» *1/2 Half Size*

1/2 HALF SIZE



W: 280 mm



L: 310 mm



Lid (Perforated)



Box (Non Perforated)

h:110 mm	 K6-500E	 K6-501E	 K6-502E	 K6-503E	 K6-504E
h:130 mm	 K6-506E	 K6-507E	 K6-508E	 K6-509E	 K6-510E
h:150 mm	 K6-512E	 K6-513E	 K6-514E	 K6-515E	 K6-516E
h:200 mm	 K6-518E	 K6-519E	 K6-520E	 K6-521E	 K6-522E



STANDARD MODEL
29

BIO-BARRIER MODEL
63

Polymer's PPSU Lids
97

WINDOW MODEL
119

B₃ Barrier MODEL
139

PLASMA MODEL
173

A1 MODEL
185

Dental, Mini, Special
205

BASKETS
223


Autolock Mechanism Emulating Indicator
243

ACCESSORIES
251


BOX (non perforated)




outside dimension



outside dimension



inside dimension



outside dimension	inside dimension	weight
310 x 272 h:100 mm	253 x 258 h:94 mm	1267 g
310 x 272 h:116 mm	253 x 258 h:110 mm	1315 g
310 x 272 h:138 mm	253 x 258 h:132 mm	1427 g
310 x 272 h:190 mm	253 x 258 h:184 mm	1660 g

PERFORATED LID



Dimension: 312 x 280 mm
Weight : 722 g

BASKET MADE FROM STAINLESS STEEL SHEET

BASKET (with stand)



	code	outside dimension
	K3-752	244 x 253 h:40 mm
	K3-753	244 x 253 h:60 mm
	K3-754	244 x 253 h:80 mm
	K3-755	244 x 253 h:110 mm

➤ The baskets which are without lid and less than 44 mm in height do not include the handles while the ones with lid include the handles.

STAND



code	height
K3-767	20 mm



STAND



code	height
K3-768	10 mm



STAND



code	height
K3-766	5 mm



BASKET MADE FROM STAINLESS STEEL WIRE MESH




code	outside dimension
K2-490	246 x 246 x h:30 mm
K2-491	246 x 246 x h:50 mm
K2-492	246 x 246 x h:70 mm
K2-493	246 x 246 x h:100 mm

SILICONE MAT



K2-521
250 x 235 mm

SILICONE MESH MAT



K2-530
250 x 240 mm

SILICONE HANDLE



K2-710	
K2-711	
K2-712	
K2-713	
K2-714	
K2-715	

PAPER FILTER



K2-740
Single use 100 pcs
with indicator. Aygun mark.
• Diameter: 184 mm

ALUMINIUM FILTER RETAINER




K2-722
Diameter: 186 mm

FILTER DISPENSER




K2-730

IDENTIFICATION LABELS (aluminium)



K2-700	
K2-701	
K2-702	
K2-703	
K2-704	
K2-705	


LABELS (paper)



K2-673
100 pcs. With indicator.

For small size and
special containers

LABELS (paper)



K2-672
100 pcs. With indicator.

For full and
Half size containers

SECURITY SEAL (blue)



K2-663
100 pcs. with indicator.
made of plastic

SECURITY SEAL (blue)



K2-664
100 pcs. made of plastic

SECURITY SEAL (blue)



K2-661
100 pcs. with number.
made of plastic

SECURITY SEAL (red)



K2-662
100 pcs. with indicator.
made of plastic

SECURITY SEAL (red)



K2-660
100 pcs. made of plastic

SECURITY SEAL (red)



K2-665
100 pcs. with number.
made of plastic

STANDARD MODEL	29
BIO-BARRIER MODEL	63
Polymers PPSU Lids	97
WINDOW MODEL	119
B ₃ Barrier MODEL	139
PLASMA MODEL	173
A1 MODEL	185
Dental, Mini, Special	205
BASKETS	223
Autolock Mechanism Emulating Indicator	243
ACCESSORIES	251

A1 MODEL

1/2 half size

