

STERILE CONTAINER SYSTEMS STANDARD MODEL



1/1 FULL SIZE



3/4 QUARTER SIZE



1/2 HALF SIZE



EXTRA LARGE SIZE



EXTRA LONG SIZE

Aygün standard model container systems are designed for perfect steam flow and penetration for loads suitable to the standards and can be used with single use paper filter or reusable PTFE Teflon filters. The standard model sterilization containers are equipped with silicone coated carrying handles which have 90° mobility and its ergonomic design provides easy and secure carrying facilities. The gaskets inside the lid and filter retainer ensure perfect tightness while locking mechanism of this system provides excellent grip with its flexible structure.

Aygün standard model sterilization container systems are manufactured as five different length, width and height and can be supplied with perforated or non-perforated boxes, perforated lids and safety lids for effective external protection. These safety lids do not interfere with sterilization and they provides protection to the perforated lid and the filter system against hard objects, liquids and dust as a second lid. All models and sizes have specific baskets and accessories which can be viewed in following parts. Modular auto lock mechanism (thermo lock) promising sterility assurance and indicator window providing visibility of the indicator card result by visual control without opening the lid are also available for this sterilization container model.

STERILE CONTAINER SYSTEMS





Standard Model

» *1/1 Full Size*

1/1 STANDARD MODEL

W: 285 mm

L: 600 mm




Specifications		Lid (Perforated)		Safety Lid	
		Box (Non Perforated)	Box (Perforated)	Box (Non Perforated)	Box (Perforated)
Height	H: 110 mm				
		K6-010	K6-040	K6-150	K6-180
		K6-011	K6-041	K6-151	K6-181
		K6-012	K6-042	K6-152	K6-182
		K6-013	K6-043	K6-153	K6-183
H: 135 mm		K6-014	K6-044	K6-154	K6-184
		K6-015	K6-045	K6-155	K6-185
		K6-016	K6-046	K6-156	K6-186
		K6-017	K6-047	K6-157	K6-187
		K6-018	K6-048	K6-158	K6-188
H: 155 mm		K6-019	K6-049	K6-159	K6-189
		K6-020	K6-050	K6-160	K6-190
		K6-021	K6-051	K6-161	K6-191
		K6-022	K6-052	K6-162	K6-192
		K6-023	K6-053	K6-163	K6-193
H: 205 mm		K6-024	K6-054	K6-164	K6-194
		K6-025	K6-055	K6-165	K6-195
		K6-026	K6-056	K6-166	K6-196
		K6-027	K6-057	K6-167	K6-197
		K6-028	K6-058	K6-168	K6-198
H: 255 mm		K6-029	K6-059	K6-169	K6-199
		K6-030	K6-060	K6-170	K6-200
		K6-031	K6-061	K6-171	K6-201
		K6-032	K6-062	K6-172	K6-202
		K6-033	K6-063	K6-173	K6-203
H: 255 mm		K6-034	K6-064	K6-174	K6-204
		K6-035	K6-065	K6-175	K6-205
		K6-036	K6-066	K6-176	K6-206
		K6-037	K6-067	K6-177	K6-207
		K6-038	K6-068	K6-178	K6-208
		K6-039	K6-069	K6-179	K6-209

Ps: All container measurements cover the outer size of the lids.


BOX (perforated)




outside dimension



outside dimension



inside dimension




code	outside dimension	inside dimension	weight
K6-100	600 x 272 h:100 mm	543 x 258 h:83 mm	2216 g
K6-101	600 x 272 h:116 mm	543 x 258 h:99 mm	2362 g
K6-102	600 x 272 h:138 mm	543 x 258 h:121 mm	2520 g
K6-103	600 x 272 h:190 mm	543 x 258 h:173 mm	2873 g
K6-104	600 x 272 h:235 mm	543 x 258 h:218 mm	3290 g


BOX (non perforated)




outside dimension



outside dimension



inside dimension



code	outside dimension	inside dimension	weight
K6-105	600 x 272 h:100 mm	543 x 258 h:94 mm	1900 g
K6-106	600 x 272 h:116 mm	543 x 258 h:110 mm	2096 g
K6-107	600 x 272 h:138 mm	543 x 258 h:132 mm	2205 g
K6-108	600 x 272 h:190 mm	543 x 258 h:184 mm	2648 g
K6-109	600 x 272 h:235 mm	543 x 258 h:229 mm	3040 g

PERFORATED LIDS



code	color	weight
K6-120		1310 g
K6-121		
K6-122		
K6-123		
K6-124		
K6-125		

Dimension: 602 x 285 mm

REMOVAL LIDS (non perforated)



code	color	weight
K6-130		1046 g
K6-131		
K6-132		
K6-133		
K6-134		
K6-135		

Dimension: 602 x 285 mm

SAFETY LIDS



code	color	weight
K6-230		988 g
K6-231		
K6-232		
K6-233		
K6-234		
K6-235		

Dimension: 606 x 288 mm

STANDARD MODEL

29

BIO-BARRIER MODEL

63

Polymer's PPSU Lids

97

WINDOW MODEL

119

B₃ Barrier MODEL

139

PLASMA MODEL

173

A1 MODEL

185

Dental, Mini, Special

205

BASKETS

223

Autolock Mechanism Emulating Indicator


243

ACCESSORIES

251

BASKET MADE FROM STAINLESS STEEL WIRE MESH


WIRE BASKET



code	outside dimension
K2-460	540 x 245 h:30 mm
K2-461	540 x 245 h:50 mm
K2-462	540 x 245 h:70 mm
K2-463	540 x 245 h:100 mm
K2-450	480 x 245 h:30 mm
K2-451	480 x 245 h:50 mm
K2-452	480 x 245 h:70 mm
K2-453	480 x 245 h:100 mm



WIRE BASKET WITH LID



code	outside dimension
K3-1010	540 x 245 h:35 mm
K3-1011	540 x 245 h:55 mm
K3-1012	540 x 245 h:75 mm
K3-1013	540 x 245 h:105 mm
K3-1014	ONLY LID, for baskets
K3-1020	480 x 245 h:35 mm
K3-1021	480 x 245 h:55 mm
K3-1022	480 x 245 h:75 mm
K3-1023	480 x 245 h:105 mm
K3-1024	ONLY LID, for baskets




WIRE BASKET WITH PERFORATED METAL PLATE

WIRE BASKET



code	outside dimension
K2-461 P	540 x 245 h:50 mm
K2-462 P	540 x 245 h:70 mm
K2-463 P	540 x 245 h:100 mm
K2-451 P	480 x 245 h:50 mm
K2-452 P	480 x 245 h:70 mm
K2-453 P	480 x 245 h:100 mm



WIRE BASKET WITH LID



code	outside dimension
K3-1011 P	540 x 245 h:55 mm
K3-1012 P	540 x 245 h:75 mm
K3-1013 P	540 x 245 h:105 mm
K3-1014 P	ONLY LID, for baskets
K3-1021 P	480 x 245 h:55 mm
K3-1022 P	480 x 245 h:75 mm
K3-1023 P	480 x 245 h:105 mm
K3-1024 P	ONLY LID, for baskets




BASKET MADE FROM STAINLESS STEEL SHEET

BASKET (with stand)



	code	outside dimension
Box and Lid	K3-720	546 x 259 h:44 mm
	K3-721	546 x 259 h:64 mm
	K3-722	546 x 259 h:84 mm
	K3-723	546 x 259 h:114 mm
Box	K3-740	540 x 253 h:40 mm
	K3-741	540 x 253 h:60 mm
	K3-742	540 x 253 h:80 mm
	K3-743	540 x 253 h:110 mm
Lid	K3-760	ONLY LID, for baskets



BASKET (with stand)



	code	outside dimension
Box and Lid	K3-724	491 x 259 h:44 mm
	K3-725	491 x 259 h:64 mm
	K3-726	491 x 259 h:84 mm
	K3-727	491 x 259 h:114 mm
Box	K3-744	485 x 253 h:40 mm
	K3-745	485 x 253 h:60 mm
	K3-746	485 x 253 h:80 mm
	K3-747	485 x 253 h:110 mm
Lid	K3-761	ONLY LID, for baskets

➤ The baskets which are without lid and less than 44 mm in height do not include the handles while the ones with lid include the handles.

STAND



code	height
K3-767	20 mm

STAND



code	height
K3-768	10 mm

STAND



code	height
K3-766	5 mm

IDENTIFICATION (aluminium)

LABEL HOLDER



code	color
K2-560	Red
K2-561	Blue
K2-562	Green
K2-563	Yellow
K2-564	Grey
K2-565	Black

LABELS



code	color
K2-570	Red
K2-571	Blue
K2-572	Green
K2-573	Yellow
K2-574	Grey
K2-575	Black

Coloured identification labels
Engraved to your specification



➤ Grey colored recommended for use in Sterrad.
Color fading By effect of Hydrogen Peroxide Gas Plasma on other colors.

STANDARD MODEL

29

BIO-BARRIER MODEL

63

Polymer's PPSU Lids

97

WINDOW MODEL

119

B₃ Barrier MODEL

139

PLASMA MODEL

173

A1 MODEL

185

Dental, Mini, Special

205

BASKETS

223

Autolock Mechanism Emulating Indicator

243

ACCESSORIES

251

DIVIDERS



K2-428
540 x 250 x 130 mm
K2-429
540 x 250 x 180 mm

WIRE BASE




K2-440
545 x 255 x 30 mm

DIVIDERS




K2-425
540 x 250 x 130 mm
K2-426
540 x 250 x 180 mm

RACK (for endoscopy)



K3-1090
505 x 250 x 125 mm
• Holds 10 endoscopic instruments
• Fits into a 1/1 size 20 cm deep container
• Wire basket 1/1 50 mm
• Container 1/1 200 mm
• Silicone mat 520 x 230 mm

DRAPE RETAINER



K2-420
545 x 255 x 30 mm

DRAPE RETAINER



K2-400
545 x 255 mm

ALUMINIUM BASKET (with lid)



K3-1119
245 x 130 x 55 mm

ALUMINIUM BASKET (with lid)




K3-1120
245 x 130 x 70 mm

ALUMINIUM BASKET (with lid)




K3-1121
245 x 130 x 105 mm

SILICONE MAT



K2-526
520 x 230 mm
K2-528
475 x 230 mm

SILICONE MESH MAT



K2-531
485 x 250 mm
K2-533
534 x 247 mm

SILICONE HANDLE



K2-710
K2-711
K2-712
K2-713
K2-714
K2-715

PAPER FILTER




K2-740
Single use 100 pcs with indicator. Aygun mark.
• Diameter: 184 mm

TEFLON FILTER



K2-742
Number of sterilization cycle > 2000
• Diameter: 194 mm

FILTER DISPENSER




K2-730

IDENTIFICATION LABELS (aluminium)




K2-700
K2-701
K2-702
K2-703
K2-704
K2-705

LABELS (paper)



K2-673
100 pcs. With indicator.
For small size and special containers

LABELS (paper)



K2-672
100 pcs. With indicator.
For full and Half size containers


(-)

AUTO LOCK MECHANISM



K1-900
(2 Pcs)

EMULATING INDICATOR CARD (class 6)



K2-674
Paper 250 Pcs.



SECURITY SEAL (blue)
K2-663
100 pcs. with indicator. made of plastic



SECURITY SEAL (blue)
K2-664
100 pcs. made of plastic



SECURITY SEAL (blue)
K2-661
100 pcs. with number. made of plastic



SECURITY SEAL (red)
K2-662
100 pcs. with indicator. made of plastic



SECURITY SEAL (red)
K2-660
100 pcs. made of plastic



SECURITY SEAL (red)
K2-665
100 pcs. with number. made of plastic



DISPOSAL STAMP
K2-717
40 x 35 x 9 mm
Requested number or letter must be informed during the order.



LABEL HOLDER
K2-716
50 x 35 x 13 mm
Holder for identification labels

(-)

Assembly and disassembly method of safety lid for safe model containers.

