

STERILE CONTAINER SYSTEMS A1 MODEL



1/1 FULL SIZE



3/4 QUARTER SIZE



1/2 HALF SIZE

Aygün A1 model sterilization container systems basically have similar features with standard model sterilization container systems. It has been designed by restricting some features of standard model which will not risk any of sterilization necessities. When compared to standard model sterilization container systems there is no colorful lid option in this model and coloring parts are carried out to label plates. It can be supplied with five main colors except for black. Moreover the handles of this model are not silicone coated and this model cannot be supplied with perforated boxes and more than 200 mm deep boxes. Similar to other models, Aygün A1 model sterilization container systems are manufactured as three different length, width and height. All models and sizes have specific baskets and accessories which can be viewed in following parts.

STERILE CONTAINER SYSTEMS



A1 Model

» 3/4 Quarter Size


3/4 QUARTER SIZE


















W: 285 mm

L: 475 mm

Lid (Perforated)



Box (Non Perforated)

h: 110 mm	 K6-240E	 K6-241E	 K6-242E	 K6-243E	 K6-244E
h: 130 mm	 K6-246E	 K6-247E	 K6-248E	 K6-249E	 K6-250E
h: 155 mm	 K6-252E	 K6-253E	 K6-254E	 K6-255E	 K6-256E



STANDARD MODEL
29

BIO-BARRIER MODEL
63

Polymer's PPSU Lids
97

WINDOW MODEL
119

B₃ Barrier MODEL
139

PLASMA MODEL
173

A1 MODEL
185




Dental, Mini, Special
205

BASKETS
223

Autolock Mechanism Emulating Indicator
243

ACCESSORIES
251

BOX (non perforated)



outside dimension	inside dimension	weight
475 x 272 h:100 mm	422 x 258 h:94 mm	1609 g
475 x 272 h:116 mm	422 x 258 h:110 mm	1714 g
475 x 272 h:138 mm	422 x 258 h:132 mm	1842 g

PERFORATED LID



Dimension: 477 x 285 mm
Weight : 1009 g

BASKET MADE FROM STAINLESS STEEL SHEET

BASKET (with stand)



	code	outside dimension
	K3-748	405 x 253 h:40 mm
	K3-749	405 x 253 h:60 mm
	K3-750	405 x 253 h:80 mm
	K3-751	405 x 253 h:110 mm

➤ The baskets which are without lid and less than 44 mm in height do not include the handles while the ones with lid include the handles.


STAND



code	height
K3-767	20 mm




STAND



code	height
K3-768	10 mm



STAND



code	height
K3-766	5 mm




BASKET MADE FROM STAINLESS STEEL WIRE MESH




code	outside dimension
K2-1500	405 x 245 x h:30 mm
K2-1501	405 x 245 x h:50 mm
K2-1502	405 x 245 x h:70 mm
K2-1503	405 x 245 x h:100 mm

SILICONE MAT



K2-524
380 x 230 mm
K2-525
440 x 230 mm







SILICONE MESH MAT



K2-534
395 x 245 mm

SILICONE HANDLE




K2-710 
K2-711 
K2-712 
K2-713 
K2-714 
K2-715 

PAPER FILTER




K2-740
Single use 100 pcs
with indicator. Aygun mark.
• Diameter: 184 mm

ALUMINIUM FILTER RETAINER




K2-722
Diameter: 186 mm







FILTER DISPENSER




K2-730

IDENTIFICATION LABELS (aluminium)



K2-700 
K2-701 
K2-702 
K2-703 
K2-704 
K2-705 


LABELS (paper)



K2-673
100 pcs. With indicator.

For small size and
special containers

LABELS (paper)



K2-672
100 pcs. With indicator.

For full and
Half size containers

SECURITY SEAL (blue)



K2-663
100 pcs. with indicator.
made of plastic

SECURITY SEAL (blue)



K2-664
100 pcs. made of plastic

SECURITY SEAL (blue)



K2-661
100 pcs. with number.
made of plastic

SECURITY SEAL (red)



K2-662
100 pcs. with indicator.
made of plastic

SECURITY SEAL (red)



K2-660
100 pcs. made of plastic

SECURITY SEAL (red)



K2-665
100 pcs. with number.
made of plastic

STANDARD MODEL	29
BIO-BARRIER MODEL	63
Polymers PPSU Lids	97
WINDOW MODEL	119
B ₃ Barrier MODEL	139
PLASMA MODEL	173
A1 MODEL	185
Dental, Mini, Special	205
BASKETS	223
Autolock Mechanism Emulating Indicator	243
ACCESSORIES	251

A1 MODEL

3/4 Quarter size

