

## STERILE CONTAINER SYSTEMS PLASMA MODEL



1/1 FULL SIZE



3/4 QUARTER SIZE



1/2 HALF SIZE



DENTAL MODEL



Mini MODEL



1/2 DENTAL MODEL



ENDO MODEL



TRAY LAPAROSCOPY

Aygün plasma model container systems are designed for both steam and plasma sterilization (destroying the microorganisms by free radicals). The suitable accessories like paper filter for steam or polypropylene filter for plasma must be chosen by the user according to the need. The plasma model containers are equipped with carrying handles which have 90° mobility and its ergonomic design provides easy and secure carrying facilities. The gaskets inside the lid and filter retainer ensure perfect tightness while locking mechanism of this system provides excellent grip with its flexible structure.

Aygün plasma model sterilization container systems are manufactured with different length, width and height and can be supplied with perforated bottom and lids. All models and sizes have specific baskets and accessories which can be viewed in following parts.

## STERILE CONTAINER SYSTEMS



**LOW TEMPERATURE STERILIZATION CONTAINERS**

**For STERRAD<sup>®</sup>**



\* STERRAD is a registered trademark owned by Johnson & Johnson Company.

3/4 (QUARTER SIZE)



W: 285 mm

L: 475 mm




code	outside dimension	Specifications
<b>K6-260Z</b>	475 x 285 h:110 mm	Lid (Perforated)
<b>K6-266Z</b>	475 x 285 h:130 mm	
<b>K6-272Z</b>	475 x 285 h:155 mm	

Ps: All container measurements cover the outer size of the lids.

BOX (perforated)




outside dimension



outside dimension



inside dimension



code	outside dimension	inside dimension	weight
<b>K6-300Z</b>	475 x 272 h:100 mm	422 x 258 h:83 mm	1792 g
<b>K6-301Z</b>	475 x 272 h:116 mm	422 x 258 h:99 mm	1836 g
<b>K6-302Z</b>	475 x 272 h:138 mm	422 x 258 h:121 mm	1980 g
<b>K6-303Z</b>	475 x 272 h:185 mm	422 x 258 h:168 mm	2240 g

PERFORATED LIDS



code	weight
<b>K6-340Z</b>	1014 g


Dimension: 477 x 285 mm

SECURITY SEAL




**K2-663Z**  
100 pcs. with indicator.  
for hydrogen peroxide  
plasma sterilization

POLYPROPYLENE FILTER




**K2-741Z**  
Single use 100 pcs.  
Diameter: 184 mm  
for hydrogen peroxide  
plasma sterilization

LABEL



**K2-672Z**  
100 pcs. with indicator.  
for hydrogen peroxide  
plasma sterilization

LABEL




**K2-673Z**  
100 pcs. with indicator.  
for hydrogen peroxide  
plasma sterilization

## BASKET MADE FROM STAINLESS STEEL WIRE MESH

### WIRE BASKET



code	outside dimension
<b>K2-1500</b>	405 x 245 h:30 mm
<b>K2-1501</b>	405 x 245 h:50 mm
<b>K2-1502</b>	405 x 245 h:70 mm
<b>K2-1503</b>	405 x 245 h:100 mm
<b>K3-1034</b>	ONLY LID, for baskets



## WIRE BASKET WITH PERFORATED METAL PLATE

### WIRE BASKET




code	outside dimension
<b>K2-1501 P</b>	405 x 246 h:50 mm
<b>K2-1502 P</b>	405 x 246 h:70 mm
<b>K2-1503 P</b>	405 x 246 h:100 mm
<b>K3-1034 P</b>	ONLY LID, for baskets




## BASKET MADE FROM STAINLESS STEEL SHEET

### BASKET (with stand)



	code	outside dimension
Box	<b>K3-748</b>	405 x 253 h:40 mm
	<b>K3-749</b>	405 x 253 h:60 mm
	<b>K3-750</b>	405 x 253 h:80 mm
	<b>K3-751</b>	405 x 253 h:110 mm
Lid	<b>K3-762</b>	ONLY LID, for baskets



➤ The baskets which are without lid and less than 44 mm in height do not include the handles while the ones with lid include the handles.

### STAND



code	height
<b>K3-767</b>	20 mm

### STAND



code	height
<b>K3-766</b>	5 mm



STANDARD MODEL  
29

BIO-BARRIER MODEL  
63

Polymers PPSU Lids  
97

WINDOW MODEL  
119

B<sub>3</sub> Barrier MODEL  
139

PLASMA MODEL  
173

A1 MODEL  
185

Dental, Mini, Special  
205

BASKETS  
223

Autolock Mechanism Emulating Indicator  
243

ACCESSORIES  
251