

TethysTM



Advanced Customizing Medical Diode Laser



Wuhan Gigaa Optronics Technology Co.,Ltd.

www.gigaalaser.com

Advanced Customizing Medical Diode Laser

Tethys is the latest GIGAA LASER platform which is based on **Android system** for professional customizing machine, it satisfies the different needs of multiple disciplines.

Compatible up to 4 wavelength lasers Max. power to 120W

Laser combination:

- 532nm (4W/5W/8W)
- 810nm (10W/15W/20W/30W/60W)
- 940nm (10W/15W/20W/30W/60W)
- 980nm (10W/15W/20W/30W/60W/120W)
- 1064nm (10W/30W)
- 1210nm(15W/30W)
- 1470nm (10W/15W/20W/30W/45W/120W)
- 1940nm (6W/10W/12W/15W)

The system utilizes the high-end laser technology and **Android system** to make each wavelength can work individual and also can use all wavelength mixed. Each wavelength can adjust the parameters individual (optional)

One machine for all round demand

- Dental
- ENT
- PLDD
- EVLT
- Proctology
- Gynecology
- Dermatology
- Podiatry
- Lipolysis
- HairRemoval
- Capillary
- Physical Therapy
- VET



Wuhan Gigaa Optronics Technology Co.,Ltd has enjoyed a good reputation as medical laser device manufacturer over the world since our establishment in 2005. Our products are exported to over 87 countries now.



ADD: B8-A5, Hi-Tech Medical Device Industrial Park, #818 Gaoxin Avenue,
Wuhan 430206, China
Ph: 86-27-67848871 67848872 Fax: 86-27-67848873
<http://www.gigaalaser.com> Email: info@gigaalaser.com

