

STERILE CONTAINER SYSTEMS STANDARD MODEL



1/1 FULL SIZE



3/4 QUARTER SIZE



1/2 HALF SIZE



EXTRA LARGE SIZE



EXTRA LONG SIZE

Aygün standard model container systems are designed for perfect steam flow and penetration for loads suitable to the standards and can be used with single use paper filter or reusable PTFE Teflon filters. The standard model sterilization containers are equipped with silicone coated carrying handles which have 90° mobility and its ergonomic design provides easy and secure carrying facilities. The gaskets inside the lid and filter retainer ensure perfect tightness while locking mechanism of this system provides excellent grip with its flexible structure.

Aygün standard model sterilization container systems are manufactured as five different length, width and height and can be supplied with perforated or non-perforated boxes, perforated lids and safety lids for effective external protection. These safety lids do not interfere with sterilization and they provides protection to the perforated lid and the filter system against hard objects, liquids and dust as a second lid. All models and sizes have specific baskets and accessories which can be viewed in following parts. Modular auto lock mechanism (thermo lock) promising sterility assurance and indicator window providing visibility of the indicator card result by visual control without opening the lid are also available for this sterilization container model.

STERILE CONTAINER SYSTEMS



Standard Model


» *Extra Large Size*





















































EXTRA LARGE SIZE



W: 372 mm

L: 600 mm




Specifications ➔		Lid (Perforated)	Lid (Perforated)
Height			
		Box (Non Perforated)	Box (Perforated)
H:120 mm		K6-136  	K6-210  
		K6-137  	K6-211  
		K6-138  	K6-212  
		K6-139  	K6-213  
		K6-140  	K6-214  
		K6-141  	K6-215  
H:185 mm		K6-142  	K6-216  
		K6-143  	K6-217  
		K6-144  	K6-218  
		K6-145  	K6-219  
		K6-146  	K6-220  
		K6-147  	K6-221  

➔ Ps: All container measurements cover the outer size of the lids.

BOX (perforated)

outside dimension



outside dimension




inside dimension




code	outside dimension	inside dimension	weight
K6-400	600 x 365 h:105 mm	543 x 350 h:88 mm	2480 g
K6-401	600 x 365 h:167 mm	543 x 350 h:150 mm	3122 g

BOX (non perforated)


outside dimension



outside dimension



inside dimension



code	outside dimension	inside dimension	weight
K6-402	600 x 365 h:105 mm	543 x 350 h:99 mm	2271 g
K6-403	600 x 365 h:167 mm	543 x 350 h:161 mm	2856 g

PERFORATED LIDS
 Dimension: 602 x 370 mm
 Weight: 1076 g

K6-326
 Grey

K6-327
 Yellow

K6-328
 Green

K6-329
 Blue

K6-330
 Red

K6-331
 Black
REMOVAL LIDS (non perforated)
 Dimension: 602 x 370 mm
 Weight: 1405 g

K6-332
 Grey

K6-333
 Yellow

K6-334
 Green

K6-335
 Blue

K6-336
 Red

K6-337
 Black
STANDARD
MODEL

29

BIO-BARRIER
MODEL

63

Polymers PPSU
Lids

97

WINDOW
MODEL

119

B₃ Barrier
MODEL

139

PLASMA
MODEL

173

A1
MODEL

185

Dental, Mini,
Special

205

BASKETS

223

Autolock Mechanism
Emulating Indicator

243

ACCESSORIES


251

BASKET MADE FROM STAINLESS STEEL WIRE MESH

WIRE BASKET




code	outside dimension
K2-1701	533 x 337 h:50 mm
K2-1702	533 x 337 h:70 mm
K2-1703	533 x 337 h:100 mm



WIRE BASKET WITH LID



code	outside dimension
K3-1704	533 x 337 h:55 mm
K3-1705	533 x 337 h:75 mm
K3-1706	533 x 337 h:105 mm
K3-1707	ONLY LID, for baskets



IDENTIFICATION (aluminium)

LABEL HOLDER



code	color
K2-560	Red
K2-561	Blue
K2-562	Green
K2-563	Yellow
K2-564	Grey
K2-565	Black

LABELS



code	color
K2-570	Red
K2-571	Blue
K2-572	Green
K2-573	Yellow
K2-574	Grey
K2-575	Black


Coloured identification labels
Engraved to your specification





➡ Grey coloured recommended for use in Sterrad.
Color fading By effect of Hydrogen Peroxide Gas Plasma on other colors.

BASKET MADE FROM STAINLESS STEEL SHEET

BASKET (with stand)



	code	outside dimension
Box and Lid	K3-701	539 x 346 h:44 mm
	K3-702	539 x 346 h:64 mm
	K3-703	539 x 346 h:84 mm
	K3-704	539 x 346 h:114 mm
Box	K3-705	533 x 340 h:40 mm
	K3-706	533 x 340 h:60 mm
	K3-707	533 x 340 h:80 mm
	K3-708	533 x 340 h:110 mm
Lid	K3-709	ONLY LID, for baskets




➡ The baskets which are without lid and less than 44 mm in height do not include the handles while the ones with lid include the handles.

STAND




code	height
K3-767	20 mm

STAND



code	height
K3-768	10 mm

STAND



code	height
K3-766	5 mm

DIVIDERS



K2-1731
533 x 337 x 70 mm
K2-1732
533 x 337 x 120 mm


(-)

DIVIDERS



K2-1721
533 x 337 x 70 mm
K2-1722
533 x 337 x 120 mm

DRAPE RETAINER




K2-1712
533 x 337 x 30 mm

DRAPE RETAINER



K2-1711
533 x 337 mm

SILICONE MESH MAT



K2-535
527 x 335 mm

SILICONE HANDLE



K2-710
K2-711
K2-712
K2-713
K2-714
K2-715

ALUMINIUM BASKET (with lid)



K3-1119
245 x 130 x 55 mm

ALUMINIUM BASKET (with lid)



K3-1120
245 x 130 x 70 mm

ALUMINIUM BASKET (with lid)



K3-1121
245 x 130 x 105 mm

PAPER FILTER




K2-740
Single use 100 pcs
with indicator. Aygun mark.
• Diameter: 184 mm

TEFLON FILTER



K2-742
Number of sterilization
cycle > 2000
• Diameter: 194 mm

FILTER DISPENSER




K2-730

IDENTIFICATION LABELS (aluminium)




K2-700
K2-701
K2-702
K2-703
K2-704
K2-705

LABELS (paper)



K2-673
100 pcs. With indicator.
For small size and
special containers

LABELS (paper)



K2-672
100 pcs. With indicator.
For full and
Half size containers


(-)

AUTO LOCK MECHANISM



K1-900
(2 Pcs)

EMULATING INDICATOR CARD (class 6)



K2-674
Paper 250 Pcs.



SECURITY SEAL (blue)
K2-663
100 pcs. with indicator.
made of plastic



SECURITY SEAL (blue)
K2-664
100 pcs. made of plastic



SECURITY SEAL (blue)
K2-661
100 pcs. with number.
made of plastic



SECURITY SEAL (red)
K2-662
100 pcs. with indicator.
made of plastic



SECURITY SEAL (red)
K2-660
100 pcs. made of plastic



SECURITY SEAL (red)
K2-665
100 pcs. with number.
made of plastic



DISPOSAL STAMP
K2-717
40 x 35 x 9 mm
Requested number or
letter must be informed
during the order.



LABEL HOLDER
K2-716
50 x 35 x 13 mm
Holder for identification
labels

(-)

Assembly and disassembly method of safety lid for safe model containers.

