

STERILE CONTAINER SYSTEMS STANDARD MODEL



1/1 FULL SIZE



3/4 QUARTER SIZE



1/2 HALF SIZE



EXTRA LARGE SIZE



EXTRA LONG SIZE

Aygün standard model container systems are designed for perfect steam flow and penetration for loads suitable to the standards and can be used with single use paper filter or reusable PTFE Teflon filters. The standard model sterilization containers are equipped with silicone coated carrying handles which have 90° mobility and its ergonomic design provides easy and secure carrying facilities. The gaskets inside the lid and filter retainer ensure perfect tightness while locking mechanism of this system provides excellent grip with its flexible structure.

Aygün standard model sterilization container systems are manufactured as five different length, width and height and can be supplied with perforated or non-perforated boxes, perforated lids and safety lids for effective external protection. These safety lids do not interfere with sterilization and they provides protection to the perforated lid and the filter system against hard objects, liquids and dust as a second lid. All models and sizes have specific baskets and accessories which can be viewed in following parts. Modular auto lock mechanism (thermo lock) promising sterility assurance and indicator window providing visibility of the indicator card result by visual control without opening the lid are also available for this sterilization container model.

STERILE CONTAINER SYSTEMS



Standard Model



» *Extra Long Size*

EXTRA LONG SIZE



W: 285 mm




L: 725 mm



Specifications		Lid (Perforated)	Lid (Perforated)
		Box (Non Perforated)	Box (Perforated)
Height H: 140 mm		K6-2010	K6-2040
		K6-2011	K6-2041
		K6-2012	K6-2042
		K6-2013	K6-2043
		K6-2014	K6-2044
		K6-2015	K6-2045
Height H: 205 mm		K6-2016	K6-2046
		K6-2017	K6-2047
		K6-2018	K6-2048
		K6-2019	K6-2049
		K6-2020	K6-2050
		K6-2021	K6-2051




Ps: All container measurements cover the outer size of the lids.

BOX (perforated)



code	outside dimension	inside dimension	weight
K6-2000	725 x 272 h:127 mm	687 x 258 h:110 mm	2646 g
K6-2001	725 x 272 h:192 mm	687 x 258 h:175 mm	3136 g

BOX (non perforated)



code	outside dimension	inside dimension	weight
K6-2005	725 x 272 h:127 mm	687 x 258 h:121 mm	2446 g
K6-2006	725 x 272 h:192 mm	687 x 258 h:186 mm	2974 g

PERFORATED LIDS



Dimension: 726 x 285 mm
Weight: 1348 g

- K6-2120**
Grey
- K6-2121**
Yellow
- K6-2122**
Green
- K6-2123**
Blue
- K6-2124**
Red
- K6-2125**
Black

REMOVAL LIDS (non perforated)



Dimension: 726 x 285 mm
Weight: 1287 g

- K6-2130**
Grey
- K6-2131**
Yellow
- K6-2132**
Green
- K6-2133**
Blue
- K6-2134**
Red
- K6-2135**
Black

STANDARD MODEL

29

BIO-BARRIER MODEL

63

Polymer's PPSU Lids

97

WINDOW MODEL

119

B₃ Barrier MODEL

139

PLASMA MODEL

173

A1 MODEL

185

Dental, Mini, Special

205

BASKETS

223

Autolock Mechanism Emulating Indicator

243

ACCESSORIES

251

BASKET MADE FROM STAINLESS STEEL WIRE MESH

WIRE BASKET



code	outside dimension
K2-1801	665 x 250 h:50 mm
K2-1802	665 x 250 h:70 mm
K2-1803	665 x 250 h:100 mm



WIRE BASKET WITH LID



code	outside dimension
K3-1804	665 x 250 h:55 mm
K3-1805	665 x 250 h:75 mm
K3-1806	665 x 250 h:105 mm
K3-1807	ONLY LID, for baskets



IDENTIFICATION (aluminium)

LABEL HOLDER



code	color
K2-560	Red
K2-561	Blue
K2-562	Green
K2-563	Yellow
K2-564	Grey
K2-565	Black

LABELS



code	color
K2-570	Red
K2-571	Blue
K2-572	Green
K2-573	Yellow
K2-574	Grey
K2-575	Black


Coloured identification labels
Engraved to your specification





➤ Grey coloured recommended for use in Sterrad.
Color fading By effect of Hydrogen Peroxide Gas Plasma on other colors.

BASKET MADE FROM STAINLESS STEEL SHEET

BASKET (with stand)




	code	outside dimension
Box and Lid	K3-715	671 x 259 h:44 mm
	K3-716	671 x 259 h:64 mm
	K3-717	671 x 259 h:84 mm
	K3-718	671 x 259 h:114 mm
Box	K3-711	665 x 253 h:40 mm
	K3-712	665 x 253 h:60 mm
	K3-713	665 x 253 h:80 mm
	K3-714	665 x 253 h:110 mm
Lid	K3-710	ONLY LID, for baskets



➤ The baskets which are without lid and less than 44 mm in height do not include the handles while the ones with lid include the handles.

STAND



code	height
K3-767	20 mm

STAND




code	height
K3-768	10 mm

STAND



code	height
K3-766	5 mm

RACK (for endoscopy)



K3-1093
550 x 250 x 130 mm

- Holds 10 endoscopic instruments
- Fits into a Extra Long size 18 cm deep container
- Wire basket Extra Long 50 mm
- Container Extra Long 180 mm
- Silicone mat 340x230 (2pcs)

RACK (for kerrison)



K3-1094
220 x 250 x 160 mm

- Holds 10 kerrison punches
- Fits into a 1/2, 1/1 20 cm and Extra Long size 18 cm deep container
- Wire basket 1/2 50 mm
- Wire basket 1/1 50 mm
- Container 1/2 200 mm
- Container 1/1 200 mm
- Wire basket Extra long 50 mm
- Container Extra long 180 mm

STANDARD MODEL

29

BIO-BARRIER MODEL

63

Polymer's PPSU Lids

97

WINDOW MODEL

119

B₃ Barrier MODEL

139

PLASMA MODEL

173

A1 MODEL

185

Dental, Mini, Special

205

BASKETS

223

Autolock Mechanism Emulating Indicator

243

ACCESSORIES

251

DIVIDERS



K2-428
540 x 250 x 130 mm
K2-429
540 x 250 x 180 mm

WIRE BASE



K2-440
545 x 255 x 30 mm

DIVIDERS



K2-425
540 x 250 x 130 mm
K2-426
540 x 250 x 180 mm

ALUMINIUM BASKET (with lid)



K3-1119
245 x 130 x 55 mm

ALUMINIUM BASKET (with lid)




K3-1120
245 x 130 x 70 mm

ALUMINIUM BASKET (with lid)




K3-1121
245 x 130 x 105 mm

SILICONE MAT



K2-522
340 x 230 mm
(2 pcs)







SILICONE MESH MAT



K2-536
658 x 247 mm

SILICONE HANDLE



K2-710 
K2-711 
K2-712 
K2-713 
K2-714 
K2-715 

PAPER FILTER




K2-740
Single use 100 pcs
with indicator. Aygun mark.
• Diameter: 184 mm

TEFLON FILTER



K2-742
Number of sterilization
cycle > 2000
• Diameter: 194 mm







FILTER DISPENSER




K2-730

IDENTIFICATION LABELS (aluminium)



K2-700 
K2-701 
K2-702 
K2-703 
K2-704 
K2-705 


LABELS (paper)



K2-673
100 pcs. With indicator.

For small size and
special containers

LABELS (paper)



K2-672
100 pcs. With indicator.

For full and
Half size containers

(-)

AUTO LOCK MECHANISM



K1-900
(2 Pcs)

EMULATING INDICATOR CARD (class 6)



K2-674
Paper 250 Pcs.



SECURITY SEAL (blue)

K2-663
100 pcs. with indicator.
made of plastic



SECURITY SEAL (blue)

K2-664
100 pcs. made of plastic



SECURITY SEAL (blue)

K2-661
100 pcs. with number.
made of plastic



SECURITY SEAL (red)

K2-662
100 pcs. with indicator.
made of plastic



SECURITY SEAL (red)

K2-660
100 pcs. made of plastic



SECURITY SEAL (red)

K2-665
100 pcs. with number.
made of plastic



DISPOSAL STAMP

K2-717
40 x 35 x 9 mm
Requested number or
letter must be informed
during the order.



LABEL HOLDER

K2-716
50 x 35 x 13 mm
Holder for identification
labels

(-)

Assembly and disassembly method of safety lid for safe model containers.

