

STERILE CONTAINER SYSTEMS PLASMA MODEL



1/1 FULL SIZE



3/4 QUARTER SIZE



1/2 HALF SIZE



DENTAL MODEL



Mini MODEL



1/2 DENTAL MODEL



ENDO MODEL



TRAY LAPAROSCOPY

Aygün plasma model container systems are designed for both steam and plasma sterilization (destroying the microorganisms by free radicals). The suitable accessories like paper filter for steam or polypropylene filter for plasma must be chosen by the user according to the need. The plasma model containers are equipped with carrying handles which have 90° mobility and its ergonomic design provides easy and secure carrying facilities. The gaskets inside the lid and filter retainer ensure perfect tightness while locking mechanism of this system provides excellent grip with its flexible structure.

Aygün plasma model sterilization container systems are manufactured with different length, width and height and can be supplied with perforated bottom and lids. All models and sizes have specific baskets and accessories which can be viewed in following parts.

STERILE CONTAINER SYSTEMS



LOW TEMPERATURE STERILIZATION CONTAINERS

For STERRAD[®]

* STERRAD is a registered trademark owned by Johnson & Johnson Company.

DENTAL MODEL




code	outside dimension	Specifications
K6-1820Z	325 x 192 x 55 mm	Lid (Perforated)
K6-1826Z	325 x 192 x 70 mm	 Box (Perforated)
K6-1832Z	325 x 192 x 135 mm	

BOX (perforated)



code	outside dimension	weight
K6-1860Z	323 x 186 x 43 mm	493 g
K6-1861Z	323 x 186 x 58 mm	535 g
K6-1862Z	323 x 186 x 122 mm	745 g

PERFORATED LIDS



code	weight
K6-1880Z	638 g


Dimension: 325 x 192 mm

POLYPROPYLENE FILTER



K3-615Z
Single use 100 pcs.
Dimension: 165 x 235 mm
for hydrogen peroxide plasma sterilization

BASKET (with stand)



code	outside dimension
K3-736	275 x 175 x 53 mm
K3-737	275 x 175 x 113 mm

MINI MODEL



code	outside dimension	Specifications
K6-1920Z	315 x 142 x 45 mm	Lid (Perforated)
K6-1926Z	315 x 142 x 70 mm	 Box (Perforated)
K6-1932Z	315 x 142 x 100 mm	

BOX (perforated)



code	outside dimension	weight
K6-1960Z	315 x 132 x 32 mm	338 g
K6-1961Z	315 x 132 x 63 mm	432 g
K6-1962Z	315 x 132 x 92 mm	523 g

PERFORATED LIDS



code	weight
K6-1970Z	513 g

Dimension: 320 x 142 mm

POLYPROPYLENE FILTER



K2-761Z
Single use 100 pcs.
Dimension: 230 x 120 mm
for hydrogen peroxide plasma sterilization


BASKET (with stand)



code	outside dimension
K3-738	265 x 125 x 63 mm

1/2 DENTAL MODEL



code	outside dimension	Specifications
K6-1520Z	227 x 166 x 65 mm	Lid (Perforated)
K6-1526Z	227 x 166 x 100 mm	
K6-1532Z	227 x 166 x 130 mm	Box (Perforated)

BOX (perforated)



code	outside dimension	weight
K6-1560Z	227 x 162 x 57 mm	386 g
K6-1561Z	227 x 162 x 88 mm	463 g
K6-1562Z	227 x 162 x 122 mm	525 g

PERFORATED LIDS



code	weight
K6-1570Z	488 g

Dimension: 220 x 170 mm

POLYPROPYLENE FILTER



K2-780Z


Single use 100 pcs.
Dimension: 130 x 110 mm
for hydrogen peroxide plasma sterilization

BASKET



code	outside dimension
K3-786	191 x 130 x 34 mm

Without silicone fixation systems



code	outside dimension
K3-786K	191 x 130 x 36 mm

Basket with silicone fixation systems
Holds 12 instruments with Ø 5 mm


SECURITY SEAL



K2-663Z

100 pcs. with indicator.
for hydrogen peroxide plasma sterilization

LABEL



K2-673Z

100 pcs. with indicator.
for hydrogen peroxide plasma sterilization

STANDARD MODEL

29

BIO-BARRIER MODEL

63

Polymer's PPSU Lids

97

WINDOW MODEL

119

B₃ Barrier MODEL

139

PLASMA MODEL

173

A1 MODEL

185

Dental, Mini, Special

205

BASKETS

223

Autolock Mechanism Emulating Indicator

243

ACCESSORIES

251

ENDO MODEL



code	outside dimension	Specifications
K6-1400Z	535 x 160 x 75 mm	Lid (Perforated)
K6-1406Z	535 x 160 x 100 mm	Box (Perforated)

BOX (perforated)



code	outside dimension	weight
K6-1100Z	528 x 145 x 66 mm	676 g
K6-1101Z	528 x 145 x 88 mm	786 g

PERFORATED LIDS



code	weight
K6-1000Z	795 g

Dimension: 532 x 160 mm

POLYPROPYLENE FILTER




K2-761Z
Single use 100 pcs.
Dimension: 230 x 120 mm
for hydrogen peroxide plasma sterilization

BASKET (with stand)




code	outside dimension
K3-1054	492 x 134 x 55 mm
K3-1055	492 x 134 x 65 mm

LABEL



K2-672Z
100 pcs. with indicator.
for hydrogen peroxide plasma sterilization

LABEL



K2-673Z
100 pcs. with indicator.
for hydrogen peroxide plasma sterilization

TRAY LAPAROSCOPY



color	code	outside dimension
Grey	KP-1398Z	505 x 80 x 80 mm
	KP-1399Z	505 x 155 x 80 mm
Yellow	KP-1400Z	505 x 80 x 80 mm
	KP-1401Z	505 x 155 x 80 mm
Blue	KP-1402Z	505 x 80 x 80 mm
	KP-1403Z	505 x 155 x 80 mm
Red	KP-1404Z	505 x 80 x 80 mm
	KP-1405Z	505 x 155 x 80 mm
Black	KP-1406Z	505 x 80 x 80 mm
	KP-1407Z	505 x 155 x 80 mm
Transparent	KP-1408Z	505 x 80 x 80 mm
	KP-1409Z	505 x 155 x 80 mm